MINISTAR ®

**Bedienungsanleitung
Operation Manual
Mode d'emploi
Manual de Instrucciones
Istruzioni per c'uso**

SCHEU
Dental Technology

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This is to certify that this machine:

Machine for all applications in the dental Pressure Moulding Technique

Type:

MINISTAR

Art.-No. 3100/1 230V
Art.-No. 3101/1 115V
Art.-No. 3102/1 100V

Technical Information:

Voltage	230 V / 115 V / 100 V
Power	Max. 850 W
Working Pressure	3 bar / 43 psi
L x W x H	420 x 190 x 350 mm
Gross Weight	9 kg

is in accordance with the following guide line:

Low Voltage Guide Line	2006 / 95 / EG
EMV-Guide Line	2004 / 108 / EG

Testing according to:
DIN VDE 0411-1 / EN 61010-1 / IEC 61010-1




Christian Scheu, managing director of SCHEU-DENTAL GmbH

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

Operation Manual

- I** Place the machine on a level work surface allowing sufficient space to operate the pressure chamber ①, heater ② and locking handle ③. ⚠ Avoid placing inflammable material close to the heater.
- II** Connect the main cable ① to the corresponding network (230 V/115 V, 50/60 Hz). Two fuses are integrated under the plug ②. Connection to compressed air is done by quick connect coupling ③ to the air pressure system. If the maximum standard air pressure deviates from 3 bar, the pressure can be adjusted by the pressure regulator ④ on the back of the machine. Release the pressure regulator by pulling the knob, adjust the pressure by turning (). The air pressure can be controlled by the „bar“ button on the panel.
- III** The front panel has the following operating elements:
- ① ON/OFF switch
 - ② panel indicating heating/cooling time and air pressure
 - ③ Duo-LED: green = operating mode, red = heating mode
 - ④ „+/-“ button to input heating time in steps of 5 seconds
 - ⑤ „air“ button for depressurization
 - ⑥ blue LED for indication of cooling time (continuous light) and depressurization (flashing light)
 - ⑦ „bar“ button indicates working pressure
- The machine is activated by using the illuminated ON/OFF switch ①. The green LED indicates the machine is ready to use. Heating times can be programmed in 5 second intervals by the „+/-“ buttons. Swivelling the heater over the foil activates heating immediately (LED flashes red).
⚠ Heater and heater arm become very hot; avoid touching it.
- IV** Only round foils and plates 125 mm diameter can be placed on the adjusting ring in the pressure chamber and clamped with the bayonet of the locking ring ①. Close by turning it 50° to the left.
- V** The pellets for embedding the model are in the pellet container ① which is removable by pulling it upwards. This particular movement also pulls a metal plate ② upwards closing the deep pellet channel around the model cup.
- VI** Acoustic signals indicate the end of the heating process (+5 sec. until 0 sec. = short beeps/0 sec. until - 5 sec. = continuous signal). At 000 seconds the heater should be turned to the back, at -5 seconds the heater switches off automatically for safety reasons. Turn the pressure chamber ① 180° to the right (⚠ Make sure that the rim of the model cup is clean) and turn locking ring ② 180° to the front. The magnetic valve opens and the chamber is filled immediately. Cooling time starts automatically, which is indicated by the blue LED (continuous light).
- VII** Acoustic and optical signals indicate the end of the cooling process (blue LED flashes/4 short beeps). Press „air“ button to depressurize the chamber, 4 additional short beeps (after 3 seconds) indicate the end of depressurization. Turn locking handle ①  to the back, turn locking ring ②  to the left and open pressure chamber ③ . The machine is now ready for the next pressurizing procedure.
- VIII** **Hint:** If the chamber is closed and the locking handle opened, time intervals of 60/120/180 seconds can be programmed for fabrication of positioners, mouth protectors and splints. This is done in connection with the MINISTAR® heating support ① (REF 3451) for warming up or bonding occlusal surfaces. After the corresponding heating time is programmed, the heater with heating support is turned to the front in heating position. Depending on the thickness of the material, this can be repeated several times.

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ⓖB The cooling times can be prolonged without any problem.

If necessary, the cooling time can be interrupted by pressing continuously for 5 seconds the "air" button.

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ⓖB Due to the new short-wave infrared heater heating times may have to be reduced when pressurizing plates with black content because of higher absorption; please see hints on material labels.

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