BIOSTAR®

BIOSTAR IV 1999-2006

BIOSTAR V 2006-2009

BIOSTAR VI ab / since / à partir 2009

MINISTAR S®

MINISTAR S 2003-2011

MINISTAR Scan ab / since / à partir 2011

TWINSTAR® P

TWINSTAR P ab / since / à partir 2007

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Dear Customer,

since the foundation of our company more than 80 years ago, the high functionality of our products and their safe application have been important cornerstones of our success.

Therefore, we like to give you the following important security advice on the basis of our market investigations taking place on a regular base:

In our pressure moulding machines, pressure is built up by a magnetic valve, which is released by an electric impulse. Impurities or moisture in the magnetic valve caused by the compressed air supply, however, might prevent the magnetic valve from completely closing, so that air might be leaking into the pressure chamber. This may lead to an unwanted pressure build-up in the pressure chamber, if at the same time a pressure moulding sheet is fixed in the opened pressure chamber, which may cause a safety risk.

You can exclude any risk by following these precautions:

On the illustration of the locking ring (picture 1) you can see three flat springs. **Please remove the flat spring that is on the 11.00 h position by means of a wrench (see picture 2).** After having done so, the locking ring is only held by 2 flat springs on the pressure chamber guidance. In case of an unwanted pressure build-up in the chamber, a forced decompression is automatically released. The locking ring's function of holding the pressure moulding sheet is still completely guaranteed.

When closing the locking ring with pressure moulding sheets of more than 2 mm thickness, make sure to use both hands, effecting pressure to the centre of the chamber (see picture 3). In order to lock the pressure chamber with the shaft collar, it may be necessary to press once more against the pressure chamber (see picture 4).

Moreover, we recommend having the compressed air system in your house controlled and maintained by expert staff on a regular base, in order to avoid impurities and moisture in the system. Bearing this in mind, you can rest assured your pressure moulding machine will work properly and absolutely reliable for many years!

We thank you for your comprehension and hope you'll keep on working successfully with our pressure moulding technique. Please don't hesitate to contact our field technicians for any questions.

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