BIOSTAR®

SCHEU

Bedienungsanleitung
Operation Manual
Mode d'emploi
Manual de Instrucciones
Istruzioni per l'uso

Abb. / Fig. 1



Abb. / Fig. 2

Abb. / Fig. 3

Abb. / Fig. 4

Technische Daten / Technical data

Spannung / AC	230 V, 115 V, 100 V / 50/60 Hz
Leistung / Power	750 VA
B x H x T / W x H x D	520 x 260 x 440 mm
Gewicht / Weight	14 kg
Druck / Pressure	max. 6 bar / 87 psi

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1. Safety Instructions

Carefully read the instruction manual and these important notes before commissioning and connecting the device, and make sure that the requirements are fulfilled.

General notes

- Inspect the machine for damage before commissioning. If any damage is found, it must not be operated.
- Always supervise the machine when it is switched on.
- Only use blanks of a diameter of 125 mm!
- Keep the device out of the reach of children!
- In order to protect the device against electric shock, any contact with liquids has to be avoided.
- Infrared heater and heater arm reach high temperatures during operation.
- Avoid flammable material coming close to the infrared heater.
- Only touch the handles, to avoid any risk of burns.
- Transport of the machine has to be effected solely if the machine is cooled down.
- Make sure there the heater stays free accessible during transport to avoid any damaging.
- Safe operation of the device can only be ensured if original SCHEU-DENTAL spare parts are used. Servicing must always be performed by qualified SCHEU-DENTAL service technicians or at SCHEU-DENTAL. Conversions or modifications to the device might impair product safety and lead to loss of warranty claims.
- Never open the device without authorisation. **Caution: Risk of injury!**
- Never carry out repairs on your own!
- The device is designed for power supplies with an alternating voltage of 100 V, 115 V or 230 V. (Please observe the serial number label on the rear of the device).

Installation of power supply on site

The installation of the power supply for connection of our devices must take place in accordance with international provisions and the regulations derived thereof.

Our devices are constructed according to VDE protection class I and must be connected to a socket with protective conductor contact.

Environment

The machine is fabricated from reusable or recyclable material. In accordance with the directive 2012/19 Waste Electrical and Electronic Equipment legislation (WEEE), SCHEU-DENTAL accepts to take back the device under the condition of free return delivery. Alternatively, the machine can be sent to a designated collection point for the recycling of waste electrical and electronic equipment.

2. General commissioning

⚠ Please make sure to understand and follow the below instructions to make optimal use of the BIOSTAR® machine. When placing the BIOSTAR® in its working environment, make sure there is sufficient space for the infrared heater ❶, the pressure chamber ❷ and the locking handle ❸. The BIOSTAR® should be placed in a clean environment which is free from dust and plaster. Do not place the machine near to a steam cleaner or sandblasting unit.

⚠ Fill the stainless steel pellets provided with the machine in the model cup ❹ and the remaining quantity in the pellet receiver ❺; when using the model platform ❻ make sure its surface in the model cup is free from pellets. The bearing and top and under the model cup have to be free from pellets as well.

⚠ When embedding, only use the stainless steel pellets provided with the machine. Other material (e.g. rice or sand) might lead to serious damage of the pneumatic part and the air outlet. If lighter weight materials are used these may ingest into the air stream during pressurising and evacuation. Connect the BIOSTAR® machine by plugging in only the original power and air connections; make sure the compressed air supply used is free from oil residues and humidity and has a minimum feed pressure of 6.0 bar.

3. System settings

Customise the settings of language, pressure, temperature and screen contrast according to your individual needs.

Switching on

Press 5 times the button „CE+home“ when the display of the actual machine version appears on the display.

Language selection

Select figure – confirm pressing button „code/scan“

Pressure selection

Select figure – confirm pressing button „code/scan“

Temperature selection

Select figure – confirm pressing button „code/scan“

Contrast selection

Select figure – confirm pressing button „code/scan“

4. Info menu

The following information is given in the info menu:

Software: Display of the installed firmware version

Safety instructions: Here you'll find important instructions on the safe usage of your pressure moulding device.

Indications: Discover the variety of applications in the digitised version of the application booklet „pressure moulding technique“.

Guideline to solve technical problems: Find quick and easy solutions to problems in the working process.

Range of pressure moulding material: Get an overview in the wide range of pressure moulding material and information on each kind of material.

Videos: Update your knowledge watching our short videos on the fabrication process of various kinds of pressure moulded splints.

By pressing the button „info“ **20** the info menu can be called up. The machine has to be in the basic position (pls. refer to point 7).

You can navigate through the system pressing the following buttons:

The LEDs at the buttons indicate that information is available.

Info button 20: Calling up the info menu

Hint!

After having scanned pressure moulding material, only specific information on this kind of material will be available. The info menu is not available and disappears during the pressure moulding process (heating/pressurising/cooling). Hints and instruction on the pressure moulding process are considered with priority and appear on the display, as soon as infrared heater, pressure chamber or locking ring are moved from the basic position.

← and → buttons: The horizontal navigation keys are only available after having called up an „indication“ and enable you to navigate through the various „indications“.

OK button 22: Confirming the selection.

↑ and ↓ buttons: The vertical navigation keys will lead you through the actual menu.

↶-Button 21: By pressing the **back button**, you'll return to the previous selection.

CE+home button: By pressing the button „CE+home“, you'll return to the display: „material selection depending on indication“.

5. Main connection

Connect the BIOSTAR® machine to the corresponding network by plugging in the original power cable on the back of the machine **7**, (please also refer to connection to the local network in the safety instructions). Two fuses are built in to protect the machine from high voltage (230 V version). To exchange the fuses press the little clip below the connector **8** and take off the fuse. For the 100 V/115 V versions, these fuses are built in two separate fuse covers **18**, which can be unscrewed.

6. Compressed air connection

A 20 bar tube **9** with filter element is fixed to the back of the machine and should not be exchanged. Connect the BIOSTAR® machine to the air pressure system. Use the 20 bar tube and optionally a quick connect coupling. The working pressure is set up at 6 bar/87 psi and should not be changed as it is calibrated to achieve the best pressure moulding results.

However, if you wish to lower the working pressure, just adjust it by using the pressure regulator **10** on the back of the machine. Before turning the regulator, pull to release it (↶↷).

Adjust the pressure by turning the regulator and push it in afterwards. If the pressure exceeds 6.5 bar/94 psi the security valve will be set into action. (depressurisation noise).

7. Basic position

Infrared heater **1** is swivelled back to the stop, pressure chamber **2** is opened to the left, locking ring handle **3** is swivelled back to the stop (12 o'clock position). Please refer to fig.1

8. Switching on

Press the ON/OFF button **13**, allow the machine 10 seconds to set itself up before using. Now the display shows the standard dialogue which is the starting point for any pressure moulding operation. Programming is done through display **15** or sanner **16** and shown in the display **14**.

9. Positioning of blanks

All round blanks and foils with a diameter of 125 mm and a thickness of 0.1 to 5.0 mm can be pressed in the BIOSTAR® machine. Open the pressure chamber **2** by swivelling to the left side, take off the locking ring **11** by turning it in 5 o'clock position. Place the material directly on the O ring in the pressure chamber. Put on the adjustment ring **12**, so that the bayonet clamps slip under the bevel of the locking device of the chamber. Tighten by turning clockwise into 6 o'clock position.

10. Model preparation

When using hard or hard-elastic material, models are placed into the model cup ④, filled up to its upper rim with stainless steel pellets. When using soft elastic material, models can be placed on the working platform ⑥. In both cases, the model should be placed with the incisors facing to the right towards the locking handle ③ or to the left towards the pressure chamber ② to ensure equal thickness of the material between the quadrants.

11. Scanning/programming by setting code or time

Our pressure moulding material has a 3-digit-code where all required information on temperature, heating time and cooling time is programmed. This code is printed in plain text and as bar code on our pressure moulding blanks or foils, except for material with thicknesses below 0.1 and 0.15 mm. After selecting your material, press the button „code/scan“.

The scanner ⑬ is now activated. Present only now the bar code of the chosen material to the scanner until you hear a beep. Alternatively, the code can be programmed manually.

Program the digits and confirm by pressing the button „code/scan“ again. You either have the option to enter the heating times manually. Doing so, press the key with the clock symbol ⑰, enter the required heating time and confirm by pressing the button ⑰ again.

12. Heating

After programming, directly swivel the infrared heater ① over the fixed material. The temperature is regulated automatically by a thermo element. The colour of the heater might vary from bright to dark.

13. Pressurising

Before the end of the heating cycle is reached, audible beeping sounds are heard (+ 5 sec. up to 0 sec. = beeping sounds/ 0 sec. up to 5 sec. continuous beeping sound). The heater should be swivelled back at 0 sec., at -5 sec. the heater is switched off automatically for safety reasons. Bring the pressure chamber ② over the model by turning it 180° to the right (**make sure the sealing surface is clean**) and lock it into place by turning the locking handle ③ 180° forwards. Automatically the magnetic valve opens, the pressure chamber fills with air and the cooling time starts, signal by blue LED.

14. Depressurising

Optical signals (blue flashing LED) and acoustic signals (4 beeping sounds) indicate the cooling time has elapsed. Press the flashing button „air“ to depressurise the pressure chamber. After approx. 3 seconds 4 beeping sounds indicate the end of depressurising. Now turn the locking handle ③ 180° backwards, push the locking ring ⑪ ← to the left and open the pressure chamber ② ↺. Remove the model. The machine can now be used for the next pressure moulding operation.

15. Special programs using the heating support

If you want to fabricate positioners, mouth guards or splints, you can manually programme time intervals of 60/120/180/240 seconds. Doing so, the pressure chamber has to be closed and the locking handle opened. These intervals are chosen to heat up or to seal the occlusal surfaces using the heating support. (REF 3452). After having entered the corresponding heating time, the infra red heater with the heating support hooked in is swivelled forwards in the heating position. Depending on the material thickness, this operation can be repeated several times.

16. Service/Maintenance

The BIOSTAR® machine is maintenance free. However, the air tube including the filter element should be examined regularly. If moisture or debris is released, the filter must be replaced. The model cup ④ with air holes and model platform ⑥ can be cleaned with an air nozzle or steam jet. The reflecting parts of the pressure chamber should be cleaned regularly to obtain the best heating results. For cleaning remove the aluminium adjusting ring ⑫ from the pressure chamber ② completely. The sealing surfaces, the sealing ring and the 4 rivets and springs should be cleaned and lightly greased also. Place the adjusting ring back on the rivets and springs ensuring that the sealing ring comes to the top and that the whole unit is movable. Please also refer to the enclosed „Guideline to solve technical problems“.

17. Service updates

You'll receive a registration card with the machine. By signing-up, you'll be provided automatically with information on updates. You can download the updates on our website www.scheu-dental.com/downloads/support.

Then proceed to unpack the download file and copy the subdirectories „data“ and „fw“ on an empty and unformatted USB flash drive. If you need further help, click on the video under „support“. Connect the USB device to the **service interface** ⑱ on the back of the machine. Start the machine, then press 5 times the **button „OK“**. Now you'll be asked to enter the pass word for the service menu.

Enter the **figure „1“** and confirm by pressing the **button „OK“**. By pressing again the **figure „1“** the update starts. After the update is finished, quit the service menu by pressing the **button „CE+home“**.

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Due to a safety-related technical modification the pressure chamber of your pressure moulding machine has been equipped with a drill for decompression.

Locking the closed pressure chamber by means of the shaft collar, will tighten this decompression drilling. Please make sure the sealing surface is always lubricated, in order to ensure easy locking and less wear of the sealing.

A worn sealing can be easily replaced by pulling it out from the bush support (BIOSTAR®).

Before inserting a new sealing, please lubricate the exterior side of the sealing, too (picture 4).

Your accessories set contains 3 replacement seals BIOSTAR® 10 mm sealing, REF 3000.856 and silicone grease as well as pellet brush, REF 3005.800

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Guideline to solve technical problems

Error	Probable reason	Solution
The foils are partially burned or too soft	You chose a wrong code or heating-time for the foil.	Please consider the hints and data given on the label of each package.
Pressure moulding results are not satisfying	You chose a wrong code or heating-time for the foil.	Please consider the hints and data given on the label of each package.
No pressure build-up	<p>When embedding the model in pellets, the model cup has not been filled up to the closing area. So the foil slips away from the foil-fixation because it is extended too much. Under these circumstances a correct sealing is not given. Please make sure that only foils with a diameter of 125 mm are used!</p>	 <p>Fill the model cup with pellets up to the closing area.</p>
	<p>Balancing ring is not straight or not running well due to debris.</p> <p>Maintenance hint: If you follow the maintenance advices given on the opposite side you ensure the mobility of the balancing ring. Following the maintenance advices given on the opposite side, please make sure the balancing ring is mobile at any time.</p>	 <p>Pull out balancing ring of pressure chamber by pressing on one side and lifting up on the opposite side at one time.</p> <p>Now the 4 nuts with pressure springs can be pulled out. After cleaning and greasing with vaseline, reinsert ring.</p> 
Difficulties in closing pressure chamber	Pellets on the closing area of housing.	Before closing the pressure chamber make sure that the closing area is clean!
Difficulties in opening pressure chamber	The different steps for opening the pressure chamber have not been carried out correctly.	<p>After pushing the air button for depressurization and waiting for 5 seconds the following operations have to be done:</p> <ol style="list-style-type: none"> ① Turn the locking handle 180° to the back ② Turn the bayonett closing ring to the left ③ Open the pressure chamber 

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ⒼB Warranty terms and conditions for the following machines of SCHEU-DENTAL GmbH

The machine you bought has been carefully checked and found to be in perfect condition, meeting highest technical demands.

SCHEU-DENTAL GmbH does not provide warranties for:

- repair services effected by the buyer or unauthorised staff
- in case of damages due to inappropriate use or handling
- in case of a failure to comply with the handling, cleaning and maintaining prescriptions
- in case of inappropriate connection or incorrect installation
- in case of force majeure or any other influences beyond the control of SCHEU-DENTAL
- in case of any modifications to the machines
- in case of using spare parts of other manufacturers
- in case of signs of wear due to age, not affecting the main functions of the machines.

Utilisation of warranties does not extend or restart the standard warranty period. In case of unjustified utilisation of our warranties we reserve to charge a fee for the corresponding services.

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Managing Director of SCHEU-DENTAL GmbH