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
GB: REF 5309 / D: REF 5307 / ES: REF 5318

How to fabricate a CA[®] CLEAR ALIGNER splint using the CA[®] Distalizer Kit acc. to Dr. Echarri

This is achieved using the CA[®] telescope sleeve with the CA[®] guiding rod for buccal insertion and the VECTOR[®] 350 for palatal insertion.

Hint:

The pretension of the pressure spring can be individually adjusted. If the tension is too high, you can reduce the spring using a reamer, if the tension is too low, the spring can be slightly pulled apart.

Please also refer to the following technical literature:

CA[®] Book volume 2 Advanced (chapter 3: CA[®] Distalizer Kit + VECTOR[®] 350)

English version: REF 5309 / German version: REF 5307 / Spanish version: REF 5318

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CA® Distalizer Kit according to Dr. Echarri (REF 5319)

All components are also available separately:


- REF 5328.800 CA® pressure spring, 1 pc.
- REF 5327.800 CA® guiding rod, 1 pc.
- REF 5323.810 CA® telescope sleeve, 1 pc.
- REF 5326.800 Silicone tube, 1 pc.
- REF 2509 VECTOR® 350 (incl. key for expansion screws)

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
CA® guiding rod,
silicone tube,
CA® telescope sleeve

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CA® guiding rod and silicone tube*
have to be individually adjusted

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
Proceed with bending a loop for
retention on one side of the
CA® guiding rod.
Push the silicone tube onto the
other end and insert it into the
CA® telescope sleeve.
Make sure that approx. 1 mm
between sleeve and rod are
covered by the silicone tube.

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Apply pink wax to the contact point
between silicone tube and wire.

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Block out model using e.g. SIL-KITT
red (REF 3442).
Pressurise ISOFOLAN® foil
(REF 3207), reduce it at model
bottom and insert occlusal cuts.

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Then proceed with embedding
the model in the pellets and
pressurising the DURASOFT® pd
foil.

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Roughly reduce foil without taking
it off the model.

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
Separate the pressure moulded DURASOFT® pd splint using the manual saw (REF 5372) on the model, inserting the saw cut interdentally/mesially in front of the teeth you intend to distalise.

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Apply pink wax to the saw cut and remove wax residues from the foil.

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
Reduce the lower part of the VECTOR® 350 screw using the diamond cutting disc (REF 5288).

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Hint: Please note that you risk squeezing the space if you use scissors or a reamer to reduce it.

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Make sure the insertion angle for the CA® guiding rod is approx. 10° (see diagram), in order to prevent the potential risk of an open bite during the distalisation process. We recommend using the ECHARRI Ceph Template.

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Use blocking out putty (REF 3220) to protect the occlusal surface. Subsequently, the VECTOR® 350 screw and the assembled CA® telescope sleeve (pls. refer to fig. 3) should be aligned parallel to each other and fixed using a small amount of mixed BIOCRYL® RESIN (REF 3218/3221). (Polymerization in a pressure pot isn't necessary.)

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Apply DURASPLINT® (REF 3258/3259) to the components before the polymerisation in a pressure pot.

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After polymerization, remove blocking out putty and take off the upper part of the space holder using the CARAT® flat pliers (REF 2041) or similar pliers. Afterwards prefinish the appliance.

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Prepare the occlusal saw cut using the manual saw (REF 5372).

Attention:
Ensure not to touch metal parts!

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Remove the splint from the model. Separate it completely from bottom to metal parts using the manual saw.

Attention!
Make sure not to cut through the CA® guiding rod!

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Unscrew and dismantle the VECTOR® 350 screw completely.

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Now the silicone tube should be exposed in the vestibular area and then removed..

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
Round off cutting rims of the splint using DIMO®PRO Slim (REF 3376) at max. 10.000 rpm.

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
Push CA® pressure spring onto CA® guiding rod.

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Proceed with reassembling the VECTOR® 350 screw, making sure the CA® guiding rod slides back into the CA® telescope sleeve simultaneously. Verify the correct length of the CA® pressure spring following the hint on page 2.

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Now fixing the spring mesially in relaxed condition with a small amount of SD-Cyano Veneer Fast (REF 3701) for safety reasons. We recommend using the SD Pipette (REF 3703) for easy precision application.

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Finish and polish the CA® Distalizer splint as usual using the accessories as recommended on page 3.

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
Important:

Make sure not to exceed the maximal tooth movement in the CA® CLEAR ALIGNER technique acc. to Dr. P. Echarri (CA® Book volume 2, Advanced, pages 52/53).

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Completed CA® Distalizer splint acc. to Echarri after polishing.

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
Activation hint: VECTOR® 350

 1/4 rotation: 0.09 mm
 4/4 rotation: 0.35 mm
 Activation: approx. 3.00 mm

The CA® Distalizer splint can be used for single or double-sided molar distalisation. Make sure not to activate the screws simultaneously on the left and the right side in order to prevent the front teeth from tipping in the course of the double-sided distalisation.

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→ www.scheu-dental.com → SCHEU-ACADEMY → Themen: CA® CLEAR ALIGNER

Hint:

Please note this instruction does not replace attending our technical CA® courses.

We also recommend participating in the following additional course:

„Course No. 17 - module 5: CA® CLEAR ALIGNER splints with specific VECTOR® expansion screws“.

For further information and details on the course please visit our website:

→ www.scheu-dental.com → SCHEU-ACADEMY → Topics: CA® CLEAR ALIGNER