dura<sup>®</sup>  
splint LC

# Perfektionist transparent.

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**Instruction for use of  
DURASPLINT<sup>®</sup> LC**

The transparent light-curing resin  
specially developed for adjustment  
of DURAN<sup>®</sup> splints.

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We recommend the above mentioned items for the fabrication of adjusted occlusal splints made of DURAN® to facilitate your work and to ensure a user-friendly production process.

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#### Hints regarding safety at work

All skin contact with the uncured resins should be avoided. Make sure to wear gloves while working with the material. Make sure to wear protective clothing (work coat, protective goggles, medical face mask). Switch on suction unit.

#### Hint regarding light polymerisation

Light polymerisation requires a wavelength between 350 nm and 450 nm. The LC-6 Light Oven offers this spectral range.

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Articulate models in mean-value articulator following the dentist's instructions.

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Measure the model using the parallelometer taking into account the direction of insertion and mark the prosthetic equator.

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Mark splint extension in the area of slight undercuts in red.

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We recommend blocking out strong undercuts, like space closures or bridges, using BLUE-BLOKKER® or SIL-KITT.

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This will help to remove the finished splint.

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Proceed first with pressurising ISOFOLAN® and then DURAN® 1.5 mm or 2.0 mm following the instructions.

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**Important:** It is essential to use DURAN® as basis for the light-curing material.

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Reduce the pressure-moulded DURAN® roughly by means of the HM carbide bur or the SD-Disc.

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**Hint:** For optimum fit remove splint only after light-polymerised and shortening the DURASPLINT® LC material.

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The material surface in the complete splint extension should be sandblasted using aluminium oxide ( $Al_2O_3$ ) or roughened with a coarse disk DIMO® (brown, coarse).

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Take off grease residues from roughened area with ethyl alcohol to guarantee a good compound. Apply one thin layer of LC-Primer on the complete splint extension and proceed with curing for 5 minutes in LC-6 Light Oven.

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Apply LC-Insulation to the antagonist model with a brush, to prevent DURASPLINT® LC material from sticking to the model.

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Adapt one bar of the DURASPLINT® material on the pressure moulded splint. Additionally, CLEAR-BLOKKER® can be applied as modelling aid to reach a smooth result.

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**Important:** Avoid any skin contact with the uncured material. Make sure to wear gloves while working with the material.

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For easier modelling moisten your fingertips with LC-Insulation. Then proceed with modelling the material, giving it the desired shape.

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**Hint:** We recommend applying a thin layer of the material exceeding the splint extension, in order to prevent the DURASPLINT® LC build-up from scaling off as a result of material overlapping with DURAN®.

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Insert movement patterns and impressions using the articulator. Movements should be made only once in the corresponding direction. Then reopen the articulator, bringing it back in the centric position without contact to the material.

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**Hint:** Repeated movements or backward movements with the articulator closed might cause distortion due to the slightly sticky material.

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Remove the model with the pressure moulded splint from articulator and place it centrally in the LC-6 Light Oven. Then proceed with curing the material 2 times for 10 minutes. If the LC-6 Light Oven is used for a serial workflow, please make sure to cool down the LC-6 Light Oven after each polymerisation period by opening the drawer.

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After light-curing you can finish the splint still being on the working model.

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The slightly sticky inhibition layer can be removed by using isopropyl alcohol.

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For initial finishing we recommend using LC-Burs A, B and C.

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After grinding in sliding surfaces and teeth cusps according to the therapeutic instructions you can proceed with finishing the splint.

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Reduce material up to the splint extension carefully by means of the HM Carbide Cutter.

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Remove splint carefully from the model, taking it off simultaneously from both sides.

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**Hint:** This is essential to avoid cracks or fissures in the splint material.

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Then proceed with finishing the rims.

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
We recommend using the DIMO® disc coarse to smooth the surfaces of the rims or optionally sandpaper (both not included in DURASPLINT® LC Kit).

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The SD-Polishing Brush, SD-Felt Wheel and SD-High-Gloss Polishing Mop are best suited for polishing in combination with SD-Polishing Paste A and B.

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
Adjusted DURAN® splint with DURASPLINT® LC build-up after polishing.

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**Hint:** We don't recommend occlusal polishing, in order to maintain the full functionality of the adjusted splint surfaces.

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You can use CLEAR-BLOKKER® in a layer of 1-2 mm for any kind of completions or corrections. Make sure the inhibition layer stays on the object. If this layer has already been removed, apply LC-Primer to the corresponding area and proceed with curing for 5 minutes. Then apply CLEAR-BLOKKER® directly from the syringe using an application cannula. Proceed with light-curing for 10 minutes, e. g. in the LC-6 Light Oven.

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### Additional information:

DURASPLINT® LC and CLEAR-BLOKKER® do not contain methyl methacrylate monomer (LC-Primer contains methyl methacrylate monomer). Please use only light curing ovens that do not develop high temperature levels during operation (we recommend the LC-6 Light Oven).

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We recommend using CETRON® Cleansing Powder for care and cleansing.

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SCHEU-DENTAL GmbH