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### Metall-Kunststoff-Verbundsystem

SD-Metal Primer

REF 9038.1

#### Anwendung

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REF 9038.1

#### Application

SD-Metal Primer is used as a bonder between acrylics / composites and all types of metal.

After preparation, the surface to be faced is sandblasted with aluminium oxide (50 µm, 2 bar) and cleaned with a steam cleaner or acetone. A thin layer of Sebond™ is applied with a brush immediately after cleaning and allowed to dry for approx. 1 min. A second layer can be applied if required.

#### Storage

Sebond should be stored at 5 - 25° C. Close the bottle properly immediately after use and protect against the effects of direct light.

#### Shelf life

The maximum shelf life is printed on the label of each bottle. Do not use after the expiry date. Each bottle contains: 5 ml

#### C+B Opaquer

C+B Opaquer is applied with a stiff, short-haired brush. Incomplete opaquer coverage negatively affects the shade. At least two coats are required to ensure complete coverage. The flow of the opaquer is improved by mixing it on a mixing pad. The paste should be applied especially thinly to retention areas. Each layer is polymerised separately.

Spektra LED	Incisal/Dentina: 30 sec	Spektra fast: 3 min
	Cervical: 1 min	Spektra 2000: 3 min
		PLC Spektra: 3 min

If the opaquer layers are applied too thickly, it prevents optimum polymerisation and weakens the bond. Extending the polymerisation time does not increase the depth of polymerisation.

#### Useful tips

The surfaces of the framework should be cleaned of any grease and polish residues. The opaquer should only be applied in very thin layers.

#### Working time for the materials

Depending on the lighting conditions 1-3 min.

#### Storage

The material should be stored at 5 - 25° C. Close the syringe properly immediately after use and protect against the effects of direct light.

#### Shelf life

The maximum shelf life is printed on the label of each syringe. Do not use after the expiry date. Each syringe contains: 3 g

#### Guarantee

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered guidelines. As our products are subject to continued development, we reserve the right to modify them. Dialog™ complies with EN ISO 10477 which does not apply to materials for use in the posterior region which is exposed to masticatory loads.

#### Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

#### Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product - please kindly contact us.

#### Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

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