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


Operation Manual

- Place the TWINSTAR® V on a level work surface allowing sufficient space to operate the heater ① and vacuuming ②.
⚠ Avoid placing inflammable material close to the heater.

- The plug ③ with two integrated fuses ④ is on the back of the machine. Connect the enclosed mains cable to the corresponding power supply.
⚠ Pay attention for correct current voltage.

- The TWINSTAR® V has the following operating elements:

- ⑤ ON/OFF switch
- ⑥ Panel indicating heating/cooling time
- ⑦ Button to activate the heater
- ⑧ Button to stop the vacuum

The machine is activated by using the illuminated ON/OFF switch. The initial mode is **---**. Push the button ⑦ for activating the heater. The display flashes during heating time for appr. 1 minute. Display **0 0 0** indicates that the machine is ready to operate. If the heater is not used within 3 minutes, the heater switches off automatically for safety reasons. Then it has to be activated again.
⚠ Heater and heater arm become very hot; avoid touching them.

- Only round foils and plates 125 mm diameter in a thickness of 0,1 – 3,0 mm can be placed on the vacuum ring ② and clamped with the bayonet of the locking ring ⑨. Close by turning to the left. The balancing ring ⑩ in the vacuum ring has to be flexible to balance the different thicknesses of the foils.

- The model cup is closed by the model platform ⑩ on which models are placed when using ISOFOLAN® or soft elastic plates such as BIOPLAST® and COPYPLAST®. When working with hard elastic foils (e.g. DURAN® or IMPRELON®), the models are directly embedded in the pellets ⑪. In the latter remove the model platform.

- Swivelling the heater ① over the foil activates the controller for the recommended heating time. After the end of the heating process, the heater should be turned to the back. Turn the vacuum ring ② 180° to the right and press it on the O-ring. The display ⑥ indicates the cooling time automatically.

- After the recommended cooling time (Minimum cooling time 20 seconds) cancel the vacuum with the ⑧ button. After 3 seconds the display indicates **---**. Then turn locking ring ⑨ to the left and open the vacuum ring ②.

- Hint:** The display **E r r** signals that the heater ① and the vacuum ring ② are not in rest position. If the vacuum is not reduced within 180 seconds, the procedure starts automatically. This is indicated by **E r r** in the display.

Take care that the reflectance areas under the vacuum ring ② are clean otherwise the heating power is reduced.

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