VECTOR® 900

Gebrauchshinweise/Hints for use

Lagerichtiger Einbau in der OK-Dehnplatte
Correct insertion in the upper expansion plate

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Lower expansion plate with inserted VECTOR® 100 and integrated guiding surface for the guiding bars of VECTOR® 900

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Correct insertion (70°-60° relating to the occlusal plane) of the guiding bars using the angle tester acc. to Adamik

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Special screw for insertion in protrusive double plates for treatment of Angle class II cases (lower jaw retrusion/distal occlusion) as well as for bilateral expansion of the upper jaw in growing children.

Compact design with integrated expansion screw, incl. security expansion screw and SD- box. For lower jaw we recommend the combined use with VECTOR® 100 (REF 2492). For fixing and controlling the insertion the angle tester (REF 2332) is best suited.

W x H x L: 14,0 x 20,0 x 13,0 mm

Max. expansion:

5.0 mm ~ 0.0875 mm / ¼ U. ~ 57 activations / 4x1/4 rot. = 0.35 mm

Contraindications

The used stainless steel alloy is highly corrosion-resistant and suited for use in the mouth. In patients, hyper-sensitive to nickel, however, the minor nickel percentage (8-10 %) may cause allergic reactions. If the allergic reaction is confirmed by the corresponding allergy test, it may be necessary or discontinue treatment.

Hints for patients

We recommend our CETRON® care products for cleansing the appliances.