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How to fabricate a splint in the pressure moulding technique using VECTOR® Piston Screws

VECTOR® piston screws are best suited for single tooth movement during orthodontic treatment. The preliminary treatment using a VECTOR® 30, VECTOR® 40 or VECTOR® 50 expansion screws is recommended. Alternatively, an Approximal Enamel Reduction (AER) with CA® Stripping Tools can be carried out.

VECTOR® Piston Screw Assortment:

Illustration	Designation	Quantity	REF	Material	CE-Marking
	Piston spring screw length 4 mm	5 pcs.	2477	1.4305 (Spring: 1.4310)	yes
	Piston spring screw length 6 mm	5 pcs.	2478	1.4305 (Spring: 1.4310)	yes
	Piston spring screw length 8 mm	5 pcs.	2479	1.4305 (Spring: 1.4310)	yes
	Serrated housing nut length 3 mm	5 pcs.	2481	1.4305	yes
	Mounting screw*	5 pcs.	2480	–	–
	Reamer*	1 pcs.	2482	–	–
	Stainless steel holder	1 pcs.	5355	–	yes
	Screw driver*	1 pcs.	2483	–	–

* No medical product and therefore not suitable for use in the mouth! **Contains nickel!**

Working Material:

REF 3207 ISOFOLAN®
 REF 3415 DURAN® 1,0 mm
 REF 3218 BIOCRYL®-RESIN monomer
 REF 3221 BIOCRYL®-RESIN Polymer

Additionally recommended material:

REF 3442 SIL-KITT red
 REF 5276 BLUE-BLOKKER®
 REF 5371 Set-Up Bur
 REF 3265 LC-Bur C
 REF 3608 SD Tray Remover
 REF 3376 DIMO® PRO slim
 REF 3460 SD Foil Scissors A
 REF 3710 SD Polishing Paste A
 REF 3711 SD Polishing Paste B
 REF 3712 SD Polishing Brush
 REF 3713 SD Felt Wheel
 REF 3714 SD High Gloss Polishing Mop
 REF 2482 Reamer

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- Block out undercut areas using for example SIL-KITT (REF 3442).
- Apply BLUE-BLOKKER® (REF 5276) to the space holders for the tooth movement (12,22).

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- Choose the position for the piston
screw and mark.

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Insert the bore hole for the
mounting screw at an appropriate
angle using the Set-Up bur
(REF 5371).

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
Adjusting the mounting screws:
It might be necessary to shorten
the mounting screws. Ensure
that after shortening the screw
does not present any sharp edges
which might damage the serrated
nut housing. The nut has to be
unscrewed up to approx. 1 mm
below the screw head.

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Pressurise ISO-FOLAN® (REF 3207)
according to manufacturer's
instructions.

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
Insert the mounting screws (with nuts unscrewed) through the ISO-FOLAN® foil and place them into the bore holes. Shorten the ISO-FOLAN® foil at model bottom and cut interdentially. Embed model into pellets, lowering the model front more deeply into the pellet cup.

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Make sure pellets are filled in as shown on the photo and the mounting screws are vertically positioned. Then fix DURAN® (REF 3415) and mix BIOCRYL® Resin (REF 3218/3221) until it has a creamy consistency. Heat material and apply a small amount of resin to the metal parts.

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Apply resin generously during the last 5 seconds of the heating time and proceed with pressurising. Leave the model under pressure in the machine for at least 10 min. Ignore any acoustic signals during this period, making sure not to evacuate nor to switch off the machine.

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Hint: The material edges in the area of the screw head should be rounded off using the DIMO® PRO Slim wheel, in order to avoid removing too much resin. (Height and length of the piston screws have to be taken into account).

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Unscrew the mounting screws and trim the resin edge using for example the SD felt wheel (REF 3713).

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
If necessary, proceed with finishing the transition between resin and serrated housing nut with the reamer (REF 2482).

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Finish and polish conventionally (please refer to the recommended accessories listed on page 2).

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Screw in the VECTOR® piston screw (4 mm). (ATTENTION: Make sure not to unscrew the screw too far from the serrated housing nut. Doing so, control the splint on the original model (without bore holes). Please refer to the below safety note.)

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Important:

Make sure not to exceed the maximal tooth movement in the CA® CLEAR ALIGNER technique acc. to Dr. P. Echarri (CA® Book volume 2, Advanced, pages 52/53).

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Safety note:

The VECTOR® piston screws don't come with a limit stop. Therefore ensure not to unscrew them completely of the serrated housing nut!

The 4 mm screw allows a tooth movement of approx. 2-3 mm. Subsequently, the 6 mm screw should be inserted, in order to make sure the wearing comfort is not impaired by protruding screw ends.

This screw allows further tooth movements of 3-4 mm. Finally, you may use the 8 mm screw for a tooth movement up to 7 mm.

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→ www.scheu-dental.com → SCHEU-ACADEMY → Themen: CA® CLEAR ALIGNER

Hint:

Please note this instruction does not replace attending our technical CA® courses.

We also recommend participating in the following additional course:

„Course No. 17 - module 5: CA® CLEAR ALIGNER splints with specific VECTOR® expansion screws“.

For further information and details on the course please visit our website:

→ www.scheu-dental.com → SCHEU-ACADEMY → Topics: CA® CLEAR ALIGNER

SCHEU-DENTAL GmbH
www.scheu-dental.com

phone +49 2374 9288-0
fax +49 2374 9288-90

Am Burgberg 20
58642 Iserlohn · Germany

DE/GB 500/02/17 G REF PM0176.01