IP-120

IMPULSFIXATOR

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Operation Manual

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

Basic principles of arc welding

Following you will find the working principle of the IP-120 Spot Welder to quickly understand application without any problem and to achieve optimal results.

The metal pieces which are to be joined should have a minimal contact point only. The arc arises between the contact points from electrical discharge. Take care that the workpieces are not firmly pressed together. During releasing the impulse, the materials (even different materials) will melt. As the spot welding pulse is extremely short, no thermal effect to the surrounding area is created. The pieces are firmly joined and can be welded in most cases without the model. Using the TeKo principle, in some cases there is no need at all for subsequent soldering.

Important for optimal discharge and high impulse:

The crocodile clips have to be oxide free and firmly contact the working pieces. The distance between clips and melting point also influences the strength of the impulse.

Initial Operation of the IP-120

- Insert plug and switch on IP-120 ①. Preferably, the machine is switched to "auto" operation ⑧.
- Socket "1/2/3" ⑥ regulate impulse duration.
- Socket "a/b" ⑦ regulate impulse power. To achieve a 50 % higher impulse, they can be bypassed with the two-way adaptor (pic. 1).
- Control lamps ③ or ④ signal loading resp. standby after every soldering procedure.
- Control lamp ⑤ signals the electric circuit has not been disconnected, i.e. a crocodile clip has to be loosened.
- If automatic standby is not used, every soldering impulse can be loaded and deleted by means of the toggle switch ②.

Impulse powers

The lowest impulse power can be released by combination "a 1", a higher impulse by combination "b 3". All intermediate impulse powers can be chosen individually.

With the two-way adaptor with connection from "a+b" (pic. 1) the impulse power can be doubled (lowest impulse: "a+b 1", highest impulse "a+b 3").

For connection of pieces having no contact, work with 3 sockets (pic. 2). Therefore one red cable each with the corresponding clamp will be connected with socket "a+b" using the two-way adaptor and the black cable either with socket 1, 2 or 3. Then soldering is achieved with the help of a bent stainless steel wire.

Operating elements

- | | | |
|---|--------------------------------|--------------------------------|
| ① On/Off switch | ⑤ Signal for safety circuit | ⑨ Signal for automatic loading |
| ② Toggle switch for loading/unloading | ⑥ Jacks for impulse duration | ⑩ Jack for main switch |
| ③ Signal for loading process | ⑦ Jacks for impulse power | ⑪ Safety fuse |
| ④ Signal showing that appliance is ready to use | ⑧ Switch for automatic loading | ⑫ Jack for foot switch |

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IP-120 with foot pedal, safety glasses, 3 flexible cables, 3 pairs of crocodile clips, 1 two-way adaptor
Crocodile clip jagged, large
Crocodile clip jagged, medium
Crocodile clip flat, small
Two-way adaptor
Complete set of cables

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Technical Data

Primary voltage
Secondary voltage
Capacity
Fuse
Board fuse
Condensor Capacity
Measurements
Weight

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This is to certify that this machine

Machine for all application in the dental Pressure Moulding Technique

Typ:
IP-120
IMPULSFIXATOR

Art.-No. 6128/1 230V from No. 04537
Art.-No. 6209/1 115V from No. 01000

Technical Information:	
Voltage	230 V / 115 V
Power	12 VA
Working Pressure	
L x W x H	30 x 20 x 20 cm
Gross Weight	3,2 kg

is in accordance with the following guide line:

Low Voltage Guide Line 73 / 23 / EWG
EMV-Guide Line 89 / 336 / EWG

is in accordance with the following guide line:

Medicine Products 93 / 42 / EWG

Christian Scheu, managing director of SCHEU-DENTAL GmbH
Iserlohn, 06.12.2001

Christian Scheu

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Christian Scheu, Geschäftsführer der Firma Scheu-Dental GmbH
Scheu-Dental GmbH