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Report in Extracts About The Study Concerning Biocompatibility of Orthodontic Products

Orthodontic treatment is normally carried out with appliances which have direct contact with the mucous membrane for quite some time. Therefore on the hygienic point of view the materials are not allowed to release any substances which may cause harmful effects. These demands are usually summarized by the word „**Biocompatibility**“.

The main task for the biocompatibility of materials is that these materials and their released substances do not show any toxic effects. Furthermore the cytotoxic effect and the effect causing irritation of the mucous membrane are of a special significance as the materials are in direct contact with the mucous membrane.

The company SCHEU-DENTAL has placed 10 different foils, frequently used in the Pressure Moulding Technique, at Prof. Dr. G. Komposch's disposal (University of Heidelberg) in order to give her the permission to do various tests considering the aspect of biocompatibility. (Please see following chart).

Number	Trade Name		Chemical Characterization
1	Biocryl „C“	2,0 x 125 mm	Polymethylmethacrylate
2	Bioplast	2,0 x 125 mm	Ethylencopolymer and Vinylacetate. The graduation 18-24-33-40 correspond to the percentage of Vinylacetate.
3	Copyplast	1,0 x 125 mm	High molecular, thermoplastic Polymerisate on the base of Ethylen
4	Durasoft	1,8 x 125 mm	Thermoplastic Polyurethan. Polycarbonate on the base of Bisphenol A
5	Hardcast	0,6 x 125 mm	Polypropylene and other Polyolefine
6	Imprelon „S“	0,75 x 125 mm	Polycarbonate (Poly(bisphenol-A-carbonate))
7	Imprelon	1,0 x 125 mm	Polyvinylchloride
8	Imprelon	3,0 x 125 mm	High molecular, thermoplastic Polymerisate on the base of Styrene
9	Imprelon „C“	0,75 x 125mm	Copolyester *
10	Menzanium	0,9 mm hart	High-grade steel (17.9 % chrome)

Here is some information given concerning the procedure of the three tests:

Cytotoxicity

For the determination of the cytotoxicity, the *Agar-Overlay-Assay* test is used. The materials have been tested in cell-culture in the so-called Agar Diffusionstest with a neutral red colouration according to the regulation DIN 13930 on local, acute and non-specific toxicity.

Effects causing irritation of the mucous membrane

Here the *HET-CAM-Test* is done on a chicken egg as the structure of an egg is very similar to the one of the mucous membrane.

Toxicity concerning gene

One of the best known methods of testing gene toxic effects is the so-called *Ames-test*. Here the mutagenity was tested on the hereditary substance of bacteria.

None of the tested materials / foils have shown any positive reactions. (Please see following chart)

No.	Agar-Diffusion Test		HET-CAM Test	Ames-Test
	decolourization index / cell destruction index	Interpretation acc. DIN 13930	Effects after 5 min.	doubling of the Rev. Number
1	0,5 / 0,5	non cytotoxic	no effects	negative
2	0,25 / 0,25	non cytotoxic	no effects	negative
3	0,5 / 0,5	non cytotoxic	no effects	negative
4	0,25 / 0,25	non cytotoxic	no effects	negative
5	0,5 / 0,5	non cytotoxic	no effects	negative
6	0 / 0	non cytotoxic	no effects	negative
7	0,25 / 0,25	non cytotoxic	no effects	negative
8	0,5 / 0,5	non cytotoxic	no effects	negative
9	0 / 0,25	non cytotoxic	no effects	negative
10	0 / 0	non cytotoxic	no effects	negative

Summary

The results clearly state negative reactions, i. e. there is no mutagenous potential. SCHEU-DENTAL foils have no harmful effects on the patient!

* Manufacturers note: available under the brandname "DURAN®"