**Anweisungen für den Gebrauch der Schablonen von Dr. Echarri**  
**Instructions on how to use Dr. Echarri's templates**  
**Instrucciones del uso de las plantillas del Dr. Echarri**

# INHALT

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# INDEX

USE OF TEMPLATES IN DIAGNOSIS	6
Use of blue template in plaster cast study	7
Study of symmetries	7
Pont & Korkhaus Index	8
Measurement of Curve of Spee	8
Evaluation of molar torque	9
Measurement of tooth rotation	10
Use of red template in plaster cast study	11
Measurement of alveolar length	11
Use of blue template in cephalometry	12
Measurement of angles	12
Use of blue template in Cephalometric VTO	14
Determination of B1 point	14
Determination of A1 point	14
Use of green template in VTO (Visual Treatment Objective)	15
Lower incisor prediction	15
Upper incisor prediction	16
Use of green template in STO (Surgical Treatment Objective)	17
Upper incisor prediction	17
Use of green template to draw the teeth	18
Drawing of the teeth	18
Use of red template in occlusal VTO	19
Occlusal VTO procedure	19
USE OF RED TEMPLATE IN FIXED ORTHODONTICS	28
Use of red template for determination of the arch form and coordination of the upper and lower arch	28
Determination of the arch form and coordination of the upper and lower arch	29
Use of red template to form the horizontal compensation curve (toe-in)	30
Horizontal curve (toe-in)	30
Use of red template to form the sagittal compensation curve (tip-back)	31
Sagittal curve (tip-back)	31
Use of red template for loops forming	32
Checking out the loops	32
USE OF BLUE TEMPLATE IN EVALUATION OF THE CLEAR ALIGNER TREATMENT COURSE	33

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


Fig. 1A, 1B, 1C:

- A) The blue template is used for plaster cast study and for cephalometry.
- B) The red template is used for occlusal VTO and for arch forming.
- C) The green template is used for cephalometric VTO.

Fig. 1A, 1B, 1C:

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For tracing and analysis we recomend the use of acetate paper available at smile-dental\* company (REF 07-0045)

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### *Use of blue template in plaster cast study*

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Figs. 2A y 2B: Superimpose the vertical millimeter line of the template over the midline of the plaster cast. If there is a difference of more than 1 mm between the distances from bicuspid and molars of both sides to the midline, the patient presents transverse asymmetry. Using the same superimposition over the midline, superimpose the horizontal millimeter line over the mesial contact point of the left or right second bicuspid. If the same line also passes through the mesial contact point of the second bicuspid on the other side, it is considered symmetry, and if not, it is considered sagittal asymmetry.

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Fig. 3: In the same way, it can be measured: A) first bicuspid to first bicuspid distance; B) molar to molar distance; C) protrusion of incisors (Pont & Korkhaus Index).

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Fig. 4: The depth of Curve of Spee on both left and right side is measured, as it is shown in the figure 4, placing the a millimeter line of the blue template on the lower incisors edge and on the last erupted molar cusp. The depth of the Curve of Spee is measured with the other millimeter line of the template.

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Figs. 5A, 5B, 5C, 5D: Placing the template over the molar cusps and observing the plaster cast from the back (Figs. 5B & 5D), molar torque is checked. The figure 5B shows the lower plaster cast in which only the labial cusps make contact with the template, which means that the posterior teeth torque is negative. The figure 5D shows upper plaster cast in which only the palatine cusps make contact with the template, which means that the posterior teeth torque is positive.

Furthermore, it should be checked if the patient presents the interferences on non-working side during the laterality movements. If the patient shows the presence of interferences, there is a negative torque in molars, and the Wilson curve is very deep.

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Fig. 6: If the labial surface of the tooth is distally rotated, it is considered as distal rotation, and if it is mesially rotated, it is considered mesial rotation. Also, the rotations can be measured using the blue template.

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Fig. 7A & 7B: Red template should be used to measure the alveolar length.

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Figs. 8A & 8B: To measure the small dimension angles, the upper part of the blue template can be used, as it is usually used any other angle measurer (for example, interincisal angle).

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Figs. 9A & 9B: To measure large dimension angles, whose lines are not crossed in the cephalometric paper, for example, the measurement of the angle between the Frankfort and GoGn plane, the blue template should be used, based on Tales theory. One of the grid lines is positioned over one of the lines of the angle, and the zero of the half-circle is placed over the other line of the angle. In this way, the angle can be measured, without drawing the lines longer until they meet, and without drawing the parallels to one of the planes.

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Fig. 10: Two millimetric lines of the grid are used to determine the position of the lower incisor, marking a point 2 mm in front of the APO plane, and 2 mm above the occlusal plane.

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Fig. 11: The millimetric lines of the grid should be used to determine the VTO A1 point with 2 mm overjet, and 2 mm overbite in respect to A1.

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## Use of green template in VTO (Visual Treatment Objective)

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Fig. 12: To draw the lower incisor in VTO, which has to be at 22° from APo, passing through the corrected B1 which is 2 mm in front of the APo and above the occlusal plane, the 22° overlay of the green template. The vertical line of the 22° overlay is positioned over the APo plane, and the other line of the 22° overlay is positioned over the corrected B1 point. In this way, first the lower incisor axis is drawn, and then the lower incisor itself.

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Fig. 13A, 13B, 13C: To draw the upper incisor in VTO, the 130° overlay of the green template is used. One of its lines is positioned over the corrected lower incisor, and the other line over the new A1 point, at 2 mm from overjet and 2 mm from overbite. A new axis of the upper incisor is drawn, and then the upper incisor itself.

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## **Use of green template in STO (Surgical Treatment Objective)**

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Fig. 14: To draw the upper incisor at 28° from the Nasion vertical, and with the same AR point, the 28° overlay of the green template is used, placing the vertical line over the Nasion line, and the other line over the AR point. In this way, first the upper incisor axis is drawn in STO, and then the upper incisor itself.

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### Use of green template to draw the teeth

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Fig. 15: In order to uniform the teeth drawings in all tracings, the upper part of the green template can be used. To draw the incisors, the template should be positioned over the incisal axis, and incisal edge. To draw the canines, the canine axis and canine cusp should be used. To draw the molars, occlusal plane should be used, as well as the distal contact point of the molar and molar axis.

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Fig. 16: Using the red template, it is checked which arch line fits best to the lower plaster cast.

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Fig. 17: Over the lower occlusogram, it is also checked which arch line fits best to the occlusogram and the treatment plan. The arch line should pass over the incisors edges, if no retrusion or protrusion treatment is planned. If the protrusion is planned, the ideal arch line should be on the labial surface of the teeth, and if the retrusion is planned, then the ideal arch line should be on the lingual or palatine surface. The line will be moved towards labial or lingual for as many millimeters as it has been decided to protrude or retrude the teeth.

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Fig. 18: The ideal arch line is drawn in red pencil over the new sheet of acetate paper, using the red template.

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Fig. 19: As the position of posterior teeth will not be changed, they can also be drawn in the same position, and the midline is marked on the plaster cast.

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Fig. 20: One after another, the teeth are traced starting from the midline, but correcting the crowding or rotations. The central incisor is also traced.

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Fig. 22: The canine is traced. If the protrusion and retrusion calculations are correct, the canine should establish correct contact point with the first bicuspid.

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Fig. 23: The other central incisor is traced.

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Fig. 24: The other lateral incisor is traced.

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Fig. 27: The lower occlusogram and the lower occlusal VTO are super imposed over the silhouette of the posterior teeth to check the changes carried out over the anterior teeth. In this way, the treatment plan can be checked, and the number of steps necessary to carry out the treatment can be calculated.

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Fig. 28: The same procedure is followed to carry out the upper occlusal VTO, also drawn in red pencil. The silhouette of the posterior upper teeth and the midline of the occlusogram are traced.

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Fig. 29: The appropriate ideal arch line is selected using the red template, taking into account the planned protrusion or retrusion.

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Fig. 30: The first part of the upper occlusal VTO is finished.

Fig. 30: Primera parte del VTO oclusal superior realizada.

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Fig. 31: The teeth are traced over the ideal arch line, starting from the midline and correcting the crowding and rotations. The central incisor is traced.

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Fig. 34: The other central incisor is traced.

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Fig. 35: The other lateral incisor is traced.

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Fig. 38: The upper occlusogram and the upper occlusal VTO are super imposed over the silhouette of the posterior teeth to check the movements that had been carried out.

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USE OF RED TEMPLATE IN FIXED ORTHODONTICS

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*Use of red template for determination of the arch form and coordination of the upper and lower arch*

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Determination of the arch form and coordination of the upper and lower arch

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Figs. 39A & 39B: The lower arch form is positioned over the lower plaster cast to see which line adapts best to the arch form. In the maxilla, the corresponding line from the upper template will be used.

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Fig. 40: The red template is used to adapt the arch.

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*Use of red template to form the horizontal compensation curve (toe-in)*

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Fig. 41A & 41B: Distocanine points are marked in the arch with a special white pen.

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Fig. 42: Dotted lines of the template are used to form the horizontal curve beginning from the distocanine point.

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Fig. 43: As a reference, when the toe-in curve is formed, the arch should finish in the middle of the distal surface of the last erupted molar.

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**Use of red template to form the sagittal compensation curve (tip-back)**

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Fig. 44A & 44B: Sagittal curve is formed according to the overlay between the upper and lower arch overlays of the red template.

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Fig. 45: When the arch is introduced into the molar tubes, the anterior part of the arch should be at the level of the incisal necks.

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### Use of red template for loops forming

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Figs. 46A, 46B, 46C: The overlays of the three most frequent CSW technique loops are drawn in the red template. After bending them, their form and size can be checked out by positioning them over the template:

- A) super imposition of the omega loop;
- B) super imposition of the "T" loop;
- C) super imposition of the helicoidal closed "L" loop.

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## USE OF BLUE TEMPLATE IN EVALUATION OF THE CLEAR ALIGNER TREATMENT COURSE

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Fig. 47: Checking the evolution: set-up and a new plaster cast.

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Fig. 48: Checking the evolution: Measuring the set-up.

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Fig. 49: Checking the evolution: Measuring the new plaster cast.

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