Indikationstabelle

	Materialempfehlung / Recommended Material		Application Chart
1. Adapterschiene	IMPRELON® „S“	3,0 mm	Adapter Splint
2. Ätzmaske	COPYPLAST®	1,0 mm	Etching Mask
3. Aufbisschiene, adjustiert	DURAN®	1,5 mm	Adjusted Splint
4. Bissnahme- und Anprobeplatte	IMPRELON® weiß / white	2,0 / 3,0 mm	Bite registration tray
5. Bleichschiene	BIOPLAST® bleach COPYPLAST®	1,0 mm 1,0 mm	Bleaching Tray
	ohne/mit unterschiedl. Platzhaltern / without/with different spacer		
6. Bracket-Transfermaske	COPYPLAST®	0,5 mm	Bracket Transfer Mask
7. CLEAR-ALIGNER	DURAN®	0,5 mm	CLEAR-ALIGNER
		0,625 mm	
		0,75 mm	
8. CMD-Schiene, adjustiert	IMPRELON® „S“ DURASOFT® DURAN®	1,0 / 2,0 mm 1,8 mm 1,5 mm	TMJ Splint, adjusted
9. Dublierformen	BIOPLAST® (Gips / plaster) COPYPLAST® (Kunststoff / resin)	2,0 mm 2,0 mm	Model Duplications
10. Gießmasken für provisorische Kronen und Brücken	COPYPLAST®	0,5 / 1,0 mm	Moulds for Temporary Crowns and Bridges
11. Implantatschiene	DURAN®	2,0 / 3,0 mm	Template for implants
12. Immediatprothese	BIOCRYL® „C“ rosa / rose	2,0 mm	Immediate Dentures
13. Individuelle Löffel	IMPRELON® klar / clear IMPRELON® natur / nature	2,0 / 3,0 mm 3,0 mm	Custom Impression Trays
14. Invisible Retainer	IMPRELON® „S“ DURAN® COPYPLAST® C	0,5 / 0,75 / 1,0 mm 0,75 / 1,0 mm 1,0 mm	Invisible Retainer
15. Isolierfolie	ISOFOLAN®	0,10 mm	Insulating Foil
16. Kappchen	HARDCAST® COPYPLAST®	0,4 / 0,6 / 0,8 mm 0,5 / 0,6 mm	Copings
17. Kfo-Platte	BIOCRYL® „C“ klar, farbig / clear, colored BIOCRYL® „M“	1,5–3,0 mm 2,0 mm	Orthodontic Plates
18. Kieferbruchschiene	BIOCRYL® „C“ IMPRELON® „S“	2,0 mm 2,0 mm	Jaw Fracture Splint
19. Kinnkappe	IMPRELON® weiß mit / white with BIOPLAST®	2,0 mm 2,0 mm	Chin Caps
20. Knirscherschiene	IMPRELON® „S“ DURAN® DURASOFT®	1,5 / 2,0 mm 1,5–3,0 mm 1,8 / (2,5) mm	Bruxism Splints
21. Modellkaschierung	Kaschierfolie / Coating Foil	0,15 mm	Coating Foil for models
22. Mundschutz	BIOPLAST® klar, farbig / clear, colored	1,0–5,0 mm	Mouthguard
23. Positioner	BIOPLAST®	1,0–5,0 mm	Positioner
24. Prothesenbasen	BIOCRYL® „C“ rosa / rose	2,0 mm	Denture Base
25. prov. Schienen	DURAN®	0,5 / 1,0 mm	Temporary Splints
26. Retainer	BIOCRYL® „C“ klar, farbig / clear, colored BIOCRYL® „M“	1,5–3,0 mm 2,0 mm	Retainer
27. Retentionsschiene	IMPRELON® „S“	3,0 mm	Retaining Splints
28. Schnarcherschiene (IST®, TAP®)	DURAN® DURASOFT®	2,0 mm 1,8 / 2,5 mm	Snoring Splint (IST®, TAP®)
29. Set-Up Korrekturschiene	COPYPLAST® C	1,0 mm	Set-Up Splint
30. Skinverpackung	Kaschierfolie und Blisterkarton / Coating foil and blister carton	0,15 mm	Skin Packaging
31. Stripkrone	DURAN®	0,5 / 1,0 mm	Strip Crowns
32. Unterziehfolie	Platzhalterfolie / Space Maintainer Foil	0,1 mm	Space Maintainer Foil
32. Verbandplatte	DURAN®	1,0 mm	Medical Tray

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Main connection

Connect the BIOSTAR® machine to 100/115/230-V-network by plugging in the original power cable on the back of the machine **7**, two fuses are built in to protect the machine from high voltage (230 V version). To exchange the fuses please press the little clip below the connector **8** and take out the fuse. For the 100 / 115 V version, these fuses are built in two separate fuse covers **18**.

Compressed air connection

A 20 bar tube **9** with filter element is fixed to the back of the machine and should not be exchanged. Connect the BIOSTAR® machine to the air pressure system. Use the 20 bar tube and optionally a quick connect coupling. The working pressure is set up at 6 bar/87 psi and should not be changed as it is calibrated to achieve the best pressure moulding results. The active pressure is permanently shown in the display. However, if you wish to lower the working pressure, just adjust it by using the pressure regulator **10** on the back of the machine. Before turning the regulator, pull to release it. Adjust the pressure by turning the regulator (↺↻+) and push it in afterwards. If the pressure exceeds 6,5 bar/94 psi the security valve will be set into action.

Switching on

Press the ON/OFF button **13**, allow the machine 10 seconds to set itself up before using. Now the display shows the standard dia-logue which is the starting point for any pressure moulding operation. Programming is done through display **14** and key pad **15**.

Positioning of blanks

All round blanks and foils with a diameter of 125 mm and a thickness of 0,1 to 5 mm can be pressed in the BIOSTAR®. Open the pressure chamber **2** by swivelling on the left side, take off the locking ring **11** by turning it in 5 o'clock position, place the foil or blank directly on the o-ring. Put on the locking ring, so that the bayonet clamps slip under the bevel of the locking device of the chamber. Tighten by turning clockwise into 6 o'clock position.

Model preparation

When working with hard or hard elastic foils, the models are placed into the model cup **4** which is completely filled with stainless steel pellets. When working with soft elastic foils, the models can be placed on the model platform **6**. In both cases, the model should be placed with the incisors facing towards the locking handle **3** to ensure equal thickness of the foil or blank between the quadrants.

Scanning / Programming

Our pressure moulding materials have a 3 digit code in which all necessary information on temperature, heating time and cooling time is programmed. This code is printed as digit and as barcode on each pressure moulding material except material of 0,1 or 0,15 mm thickness. Choose your material, press the button „code“. The scanner **16** is activated. Present the barcode of the chosen material to the scanner until you hear a beep. Alternatively the code can be programmed manually. Program the digits and confirm by pressing the button „code“ again. You either have the option to enter the heating times manually. Doing so, press the key with the clock symbol **17**, enter the desired heating time and confirm by pressing the key **17** again.

Heating

After scanning/programming, directly swivel the infrared heater **1** over the foil. The temperature is regulated automatically by a thermo element. The colour of the heater might vary from bright to dark.

Pressurizing

Before the end of the heating cycle is reached, an audible beeping sound is heard which lasts for five seconds. After five seconds, a continuous beeping sound is heard for a further five seconds with the demand for closing the pressure chamber **2** with the locking handle **3**. During these five seconds the infrared heater is returned to its rest positon. The pressure chamber is brought over the model and locked into place by turning the locking handle 180° degrees. Automatically the magnetic valve opens, the pressure chamber fills with air and the moulding process begins. You will notice that the cooling time begins as indicated on the display. You may experience some air leakage when working with thin foils, this is normal and will not affect the final result of the moulding. If the pressure moulding is not started within ten seconds, the infrared heater will switch off and the whole process has to be started again with a new foil. Remember also to reprogram the BIOSTAR® with the correct code.

Depressurizing

After cooling time has elapsed, press the button „air“ and after another 3 seconds the locking handle **3** and the locking ring **11** can be opened. After this, the pressure chamber **2** is swivelled to the left hand side, so that the model remains in the pellets or on the model platform.

Maintenance

The BIOSTAR® machine is maintenance free. However, the air tube including the filter element should be examined regularly. If moisture or debris is released, the filter must be replaced. The model cup **4** with air holes and model platform **6** can be cleaned with an air gun. The reflecting part of the pressure chamber should be cleaned regularly to obtain the best heating results. For cleaning remove the aluminium adjusting ring **12** from the pressure chamber **2** completely. This will allow you to lightly grease the sealing surfaces. The sealing ring and the four trivets and springs should be cleaned and lightly greased also. Place the adjusting ring back on the trivets and springs ensuring that the sealing ring comes to the top and that the whole unit is movable.