# IMPRIMO® LC Denture

## Produktinformation

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# IMPRIMO® LC Denture

## Product information

<b>Description:</b>	IMPRIMO® LC Denture is a light-curing resin based on methacrylate and specifically engineered for the fabrication of denture bases; photo-polymerization is done by means of image projection systems $\leq 405$ nm. IMPRIMO® LC Denture has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Denture should be used in building processes with a z-resolution of 50- 100 $\mu\text{m}$ . Printer parameters available for Asiga MAX™ and Asiga PRO 4K.
<b>Processing instructions:</b>	Make sure IMPRIMO® LC Denture is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions. Intraoral contact with uncured or partially uncured material shall be excluded. The following specifications shall be considered for the digital design of the print object: Alveolar ridge lingual $\geq 2.5$ mm (upper / lower jaw) Palatal / lingual $\geq 3.0$ mm (upper jaw) $\geq 2.5$ mm (lower jaw) Vestibular / buccal $\geq 2.5$ mm (upper / lower jaw) We recommend using the bonding agent VITA VIONIC® BOND for fixing prosthetic teeth in the denture base (for ex. teeth printed using IMPRIMO® LC Resin). It may be required to additionally roughen the surfaces. The generally accepted professional techniques and dental procedures shall be carefully followed when preparing dentures.
<b>Contra-indications:</b>	IMPRIMO® LC Denture should not be used for any purpose other than the generative CAD/CAM fabrication of denture bases. IMPRIMO® LC Denture is a medical product class IIa.
<b>Post-processing:</b>	Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device RS wash or IMPRIMO® Cure using IMPRIMO® Cleaning Liquid.
<b>Post curing:</b>	Final and complete curing of the objects is done in the curing unit RS cure or IMPRIMO® Cure. The use of a vacuum or protective gas to prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. The polymerization process in the light oven RS cure in a vacuum takes 12 minutes. Polymerisation in the IMPRIMO® Cure takes 10 minutes in a nitrogen environment. The corresponding parameters are stored in the internal programme of the RS cure unit and IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.
<b>Shelf life / storage:</b>	Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use. Any deviation from the manufacturing method described or storage conditions may affect the mechanical and optical material properties.
<b>Hazard warnings:</b>	Prolonged skin contact with uncured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Denture might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Denture. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap. <b>Attention:</b> Polymerized resins are chemically resistant, avoid stains on clothes.
<b>Disposal:</b>	According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.
<b>Hint:</b>	To be used only for the above indication and by qualified personnel. For single use. Serious incidents in connection with the use of the product must be reported immediately to the manufacturer and the competent authority!
<b>Technical data:</b>	Viscosity: approx. 0.5 Pa s Cured material (depends on irradiation unit): Flexural strength: 114 MPa Elongation at break: 8 % Solubility: 0.41 $\mu\text{g mm}^{-3}$ Water absorption: 25.8 $\mu\text{g mm}^{-3}$ Shore hardness: 85 D Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity) *These data result from measurements of a representative sample. <b>Contains:</b> functional methacrylate resins, initiators, dyes and stabilizers. <b>Material Safety Data Sheet:</b> <a href="http://www.scheu-dental.com/en/service">www.scheu-dental.com/en/service</a>

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