IMPRIMO® LC Temp It

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Processing information

Description:

IMPRIMO® LC Temp It is a light-curing resin based on methacrylate and specifically engineered for the fabrication of long-term temporaries such as crowns, bridges and mock-ups; photo-polymerization is done by means of image projection systems (385 nm). IMPRIMO® LC Temp It has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Temp It should be used in building processes with a z-resolution of 50-100 µm. Printer parameters available for Asiga MAX™ and Asiga PRO 4K.

Processing instructions:

Make sure IMPRIMO® LC Temp It is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions. Temporary restorations can be individualized using light-curing composites. The reduction of the restoration should not exceed 0.3 mm (incisal and vestibular areas). Please follow the instructions of the manufacturer of light-cured composite materials. For cementation, common temporary tooth cements (that do not contain eugenol) and adhesives can be used. The inner areas of the crown should be conditioned by sandblasting or other commonly used techniques.

Contra-indications:

IMPRIMO® LC Temp It should not be used for any purpose other than the generative CAD / CAM fabrication of long-term temporaries such as crowns, bridges and mock-ups. Bridges with a maximum of one pontic can be fabricated. When realising frontal bridges, ensure that the connection regions must have at least 12 mm² and in the case of posterior bridges at least 14 mm². The occlusal wall thicknesses must be at least 1.5 mm (central fissure) and the circumferential thickness 1 mm. IMPRIMO® LC Temp It is a medical product class IIa.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device RS wash or IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.

Post curing:

Final and complete curing of the objects is done in the curing unit RS cure or IMPRIMO® Cure. The use of a vacuum or protective gas to prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. The polymerization process in the light oven RS cure in a vacuum takes 11 minutes. Polymerisation in the light oven IMPRIMO® Cure takes 11 minutes in a nitrogen environment. The corresponding parameters are stored in the internal programme of the RS cure and IMPRIMO® Cure unit.. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.

Shelf life / storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use. Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Temp It might occur, for more details please refer to our safety data sheet on IMPRIMO® LC Temp It. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point. May cause long-lasting harmful effects to aquatic life. All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

Note:

Serious incidents in connection with the use of the product shall be reported immediately to the manufacturer and the competent authority.

Technical data:

Viscosity: approx. 3,5 Pa s

Cured material (depends on irradiation unit):

Flexural strength: 169 MPa

Elongation at break: 4 %

Solubility: 1,1 µg mm³

Water absorption: 31,1 µg mm³

Shore hardness: 80 D

Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity)

*These data result from measurements of a representative sample.

Contains: functional methacrylate resins and inorganic fillers with particle sizes from 0.4 to 3 µm.

Material Safety Data Sheet: www.scheu-dental.com/en/service

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SCHEU-DENTAL GmbH