
CLEAR ALIGNER

**SCHEU**  
SCHEU GROUP

## CA® POWER GRIP Set

Bedienungsanleitung  
Operation Manual  
Mode d'emploi  
Manual de Instrucciones

CLEAR ALIGNER

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**Pablo Echarri, D.D.S.**

The CA®POWER GRIPS were designed by Dr. Pablo Echarri, with the help of Julio Blasco drawing the first drafts and José Fernández making the first real prototypes, both of them working in Ludent, Barcelona. The next step in developing the products were made in SCHEU-DENTAL, specially by Christian Scheu and Peter Stückrad, that have worked intensively to make it happen.

The CA®POWER GRIPS were designed to improve movements with the CA®CLEAR ALIGNER technique in a more effective and rapid way, reducing the steps of treatment.

The main concept in the MFM is to apply traction forces (elastics) in extrusion and rotation treatments at home and during the nights, and combining these forces with CA®CLEAR ALIGNER & attachments treatment during the day. The helix shape design of the MDI and DMI allow the CA®CLEAR ALIGNER to control the mesio-distal inclination of the teeth offering more retention of the appliance over the teeth achieving the movements more easily.

The FALA and FALI increase the pressure of the CA®CLEAR ALIGNER in the gingival area of the teeth making possible the torque correction.

The BUB and BUM increase the vertical dimension minimizing the occlusal interferences and improving the movements of the teeth in cross-bite position.

The forms BUAL and BUAS are very useful to increase the torque during the incisors intrusion with more effective forces.

They can also maintain the intrusion forces while the patient is not using the CA®CLEAR ALIGNER, permitting molar extrusion.

All the CA®POWER GRIPS require a special protocol of lab set-up.

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**CA® POWER GRIP Set****REF 5338.1****Content: 18 silicone forms (9 versions, 2 pcs. each); including support**

CA® POWER GRIP forms enable fabrication of individual attachments for use in the CA® CLEAR ALIGNER technique.

Set contains nine different forms which can be filled with light curing composite to be directly applied on the prepared tooth.

Curing with LED lamp.

(We recommend using Light Bond - Sealant Resin [02-0033] and Light Bond Resin - Paste [02-0031] or Tetric EvoFlow®) [www.smile-dental.de](http://www.smile-dental.de)

---

**CA® POWER GRIP, 2 silicone forms each****1 MFM Multifunctional Forced Movement**

For extrusion, mesial and distal rotation for use with CA® splints and elastics.

**2 MDI Mesial / Distal Inclination**

For setting upright frontal teeth, mesial /distal inclination.

**3 DMI Distal / Mesial Inclination**

For setting upright frontal teeth; distal / mesial inclination.

**4 FALA Force Application LABIAL**

For protrusion, retrusion and torque control up to 3°. For fabrication of activation surfaces on the teeth in the gingival area with vestibular positioning - for teeth 11, 13, 21, 23, 33 , 43.

**5 FALI Force Application LINGUAL**

For protrusion, retrusion and torque control up to 3°. For fabrication of activation surfaces on the teeth in the gingival area with labial, lingual and palatinal positioning, for teeth 12, 22, 31, 32, 41, 42.

**6 BUB Build-Up Bicuspid**

Round pad for premolars (3 x 3 mm) for occlusal bite elevation (approx. 2 mm).

**7 BUM Build-Up Molar**

Oval pad for molars (3 x 5 mm) for occlusal bite elevation (approx. 2 mm).

**8 BUAS Build-Up Anterior Small**

Occlusal bite plane for intrusion of upper frontal teeth (2.5 x 2 mm)

**9 BUAL Build-Up Anterior Large**

Occlusal bite plane for intrusion of upper frontal teeth (2.5 x 4 mm)

---

**CA® POWER GRIP Support made of stainless steel, sterilisable (REF 5340.1)**

2 pcs.

**Product Description:**

The CA® POWER GRIP technique enables fabrication of individual attachments.

Nine different forms are available which can be filled with light curing composite and directly applied on the tooth which has been prepared with Sealant/Resin. Then proceed with light curing.

**Application:**

1. Preparation of teeth: Polish the tooth surface by means of a rotating brush or with polishing paste. Then proceed with rinsing and air-drying the teeth.
2. Etching: Apply etching gel or liquid on a mixed paper. Dry the teeth to be etched and apply the etching agent on the corresponding tooth surfaces, avoiding rubbing and instead dabbing it. This takes 30 seconds, in case of milk teeth 120 seconds. Avoid any contact of soft tissue and etching gel.
3. Rinsing and drying: Rinse teeth thoroughly with water spray, each tooth for 10 seconds. In case of using etching gel for 20 seconds in order to stop etching and to remove dissolved enamel particles.
4. CA® POWER GRIP Form: Place the required CA® POWER GRIP form on the supplied support.
5. Application of primer: Apply light-curing sealant/liquid on the etched enamel surface, e.g. Assure or Light Bond Sealant.
6. Application of paste: Fill the clean form completely with light-curing primer or e.g. Light Bond past medium or LCR paste. If necessary, use a small instrument in order to fill undercuts completely.  
Positioning the attachment: Position the form on the required place for the attachment by pressing it slightly.
7. Curing of paste: Curing is done within approx. 5 - 20 seconds, depending on lamp intensity and primer.
8. Removal of form: After material has completely cured, proceed with carefully removing the form moving it slightly laterally.

**Cleansing / Disinfection:**

- The used form is to be cleansed by customary mild disinfection solutions.
- Make sure the form is completely clean before disinfection. Do not use any disinfection solutions containing formaldehydes or high-concentrated bleaching agents, as these may lead to damage and discolorations. Please follow the instructions of the corresponding manufacturer.

**To be used only by expert staff!**

The above information is for the doctor's use only, not for the patient. It is the doctor who has to evaluate this information with regard to each special treatment case and to inform the patient regarding advantages, potential side effects and complications.

**HINTS:**

To be used only for the above indication and by qualified personnel. For single use. Serious incidents in connection with the use of the product must be reported immediately to the manufacturer and the competent authority!

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SCHEU-DENTAL GmbH