

MINISTAR  [®]

 SCHEU

Bedienungsanleitung
Operation Manual
Mode d'emploi
Manual de Instrucciones
Istruzioni per l'uso

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Operation Manual

Please also refer to the detailed Security Advices.

-
- I**  Place the machine on a level work surface allowing sufficient space to operate the pressure chamber ①, heater ② and locking handle ③ as well as the scanning mechanism ④. **Infrared heater and heater arm reach high temperatures.** Only touch the handles! Avoid flammable material coming close to the infrared heater.
-
- II** Connect the main cable ① to the corresponding network (230 V / 115 V / 100 V, 50/60 Hz). Two fuses ② / ②_s are integrated in the junction block to protect the machine from high voltage. Pressurised air supply to the compressor net is provided using a standard-euro-quick coupling DN 7,2. The corresponding insert sleeve is fixed at the tube's end (6 mm inside). The tube is attached by means of a hose clip to the connector block ③. If the maximum standard air pressure deviates from 4 bar, the pressure can be adjusted by the pressure regulator ④ on the back of the machine. Release the pressure regulator by pulling the knob, adjust the pressure by turning (). The air pressure can be controlled by the „air“ button on the panel.
-
- III** The front panel has the following operating elements:
- ① ON/OFF switch
 - ② panel indicating heating/cooling time and air pressure
 - ③ Duo-LED: green = operating mode, red = heating mode
 - ④ „+/-“ button to input heating time in steps of 5 seconds
 - ⑤ „air“ button for depressurization/control of working pressure
 - ⑥ blue LED for indication of cooling time (continuous light) and depressurization (flashing light)
 - ⑦ „scan“ button for activation of scanner
- The machine is activated by using the illuminated ON/OFF switch ①. When turning on the machine, make sure it is in **basic position** with the heater being swiveled backwards, the pressure chamber and locking handle opened → handle backwards. The green LED indicates the machine is ready to use. Heating times can be programmed in 5 second intervals by the „+/-“ buttons. Swivelling the heater over the foil activates heating immediately (LED flashes red). **Heater and heater arm become very hot; avoid touching it.**
-  Scanning:
For each kind of our pressure moulding material there is a special heating time recommendation, containing as well the corresponding cooling time. This recommendation is printed in plain text and as barcode on all our pressure moulding material, except for ISOFLAN®, HARDCAST®, DURASOFT® seal and space maintainer foil. For these foils, you'll find the corresponding information on the packaging labels. Choose your material, press the button „scan“. The scanner is now activated. Present the barcode of the chosen material until you hear a beep.
-
- IV** Only round foils and plates with 125 mm diameter can be placed on the adjusting ring in the pressure chamber and clamped with the bayonet of the locking ring ①. Close by turning it 50° to the left. Before fixing any material the machine has to be ready to operate and in basic position.
-
- V** The model platform can be removed by means of the pellet brush with integrated magnet.
-
- VI** The pellets for embedding the model are in the pellet container ① which is removable by pulling it upwards. This particular movement also pulls a metal plate ② upwards closing the deep pellet channel around the model cup.
-
- VII**  Acoustic signals indicate the end of the heating process (+5 sec. until 0 sec. = short beeps/0 sec. until - 5 sec. = continuous signal). At 000 seconds the heater should be turned to the back, at -5 seconds the heater switches off automatically for safety reasons. Turn the pressure chamber ① 180° to the right (**Make sure that the rim of the model cup is clean**) and turn locking ring ② 180° to the front. The magnetic valve opens and the chamber is filled immediately. Cooling time starts automatically, which is indicated by the blue LED (continuous light).
-
- VIII**  Optical (blue flashing LED) and acoustic signals (4 short beeps) indicate the end of the coolingtime. Press „air“ button to depressurize the chamber, 4 additional short beeps (after 3 seconds) indicate the end of depressurization. Turn locking handle ①  180° to the back, turn locking ring ② ← to the left and open pressure chamber ③ . After removing the model the machine is ready for the next pressure moulding procedure.
-
- IX** **Hint:** If the chamber is closed and the locking handle opened, time intervals of 60/120/180 seconds can be programmed for fabrication of positioners, mouth protectors and splints. This is done in connection with the heating support ① (REF 3452) for warming up and laminating occlusal surfaces. After programming the heating time, turn heater with heating support to the front in heating position. Depending on the thickness of the material, this can be repeated several times.

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(GB)

The cooling times can be prolonged without any problem.

If necessary, the cooling time can be interrupted by pressing continuously for 5 seconds the "air" button.

Due to the new short-wave infrared heater heating times may have to be reduced when pressurizing plates with black content because of higher absorption; please see hints on material labels.

Please refer to our "**Guideline to solve technical problems** with BIOSTAR® VI-VI, MINISTAR®, MINISTAR S®, TWINSTAR®" on our website www.scheu-dental.com/Downloads.

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Due to a safety-related technical modification the pressure chamber of your pressure moulding machine has been equipped with a drill for decompression.

Locking the closed pressure chamber by means of the shaft collar, will tighten this decompression drilling. Please make sure the sealing surface is always lubricated, in order to ensure easy locking and less wear of the sealing.

A worn sealing can be easily replaced by pulling it out from the pressure chamber (MINISTAR®).

Before inserting a new sealing, please lubricate the exterior side of the sealing, too (picture 4).

Your accessories set contains 3 replacement seals MINISTAR® 12 mm, REF 3500.856 and silicone grease as well as pellet brush, REF 3005.800

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Guideline to solve technical problems

Error	Probable reason	Solution
The foils are partially burned or too soft	You chose a wrong code or heating-time for the foil.	Please consider the hints and data given on the label of each package.
Pressure moulding results are not satisfying	You chose a wrong code or heating-time for the foil.	Please consider the hints and data given on the label of each package.
No pressure build-up	<p>When embedding the model in pellets, the model cup has not been filled up to the closing area. So the foil slips away from the foil-fixation because it is extended too much. Under these circumstances a correct sealing is not given. Please make sure that only foils with a diameter of 125 mm are used!</p>	 <p>Fill the model cup with pellets up to the closing area.</p>
	<p>Balancing ring is not straight or not running well due to debris.</p> <p>Maintenance hint: If you follow the maintenance advices given on the opposite side you ensure the mobility of the balancing ring. Following the maintenance advices given on the opposite side, please make sure the balancing ring is mobile at any time.</p>	 <p>Pull out balancing ring of pressure chamber by pressing on one side and lifting up on the opposite side at one time.</p>  <p>Now the 4 nuts with pressure springs can be pulled out. After cleaning and greasing with vaseline, reinsert ring.</p>
Difficulties in closing pressure chamber	Pellets on the closing area of housing.	Before closing the pressure chamber make sure that the closing area is clean!
Difficulties in opening pressure chamber	The different steps for opening the pressure chamber have not been carried out correctly.	<p>After pushing the air button for depressurization and waiting for 5 seconds the following operations have to be done:</p> <ol style="list-style-type: none"> ① Turn the locking handle 180° to the back ② Turn the bayonett closing ring to the left ③ Open the pressure chamber

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Ⓖ Warranty terms and conditions for the following machines of SCHEU-DENTAL GmbH

The machine you bought has been carefully checked and found to be in perfect condition, meeting highest technical demands.

SCHEU-DENTAL GmbH does not provide warranties for:

- repair services effected by the buyer or unauthorised staff
- in case of damages due to inappropriate use or handling
- in case of a failure to comply with the handling, cleaning and maintaining prescriptions
- in case of inappropriate connection or incorrect installation
- in case of force majeure or any other influences beyond the control of SCHEU-DENTAL
- in case of any modifications to the machines
- in case of using spare parts of other manufacturers
- in case of signs of wear due to age, not affecting the main functions of the machines.

Utilisation of warranties does not extend or restart the standard warranty period. In case of unjustified utilisation of our warranties we reserve to charge a fee for the corresponding services.

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GB Security Advises:

Carefully read the instruction manual and these important notes before commissioning and connecting the device, and make sure that the requirements are fulfilled.

General notes

- Inspect the machine for damage before commissioning. If any damage is found, it must not be operated.
 - Always supervise the machine when it is switched on.
 - Only use blanks of a diameter of 125 mm!
 - Keep the device out of the reach of children!
 - In order to protect the device against electric shock, any contact with liquids has to be avoided.
 - Infrared heater and heater arm reach high temperatures during operation.
 - Avoid flammable material coming close to the infrared heater.
 - Only touch the handles, to avoid any risk of burns.
 - Transport of the machine has to be effected solely if the machine is cooled down.
 - Make sure there the heater stays free accessible during transport to avoid any damaging.
 - Safe operation of the device can only be ensured if original SCHEU-DENTAL spare parts are used. Servicing must always be performed by qualified SCHEU-DENTAL service technicians or at SCHEU-DENTAL. Conversions or modifications to the device might impair product safety and lead to loss of warranty claims.
 - Never open the device without authorisation.
- Caution: Risk of injury!**
- Never carry out repairs on your own!
 - The device is designed for power supplies with an alternating voltage of 100V, 115V or 230V. (Please observe the serial number label on the rear of the device).

Installation of power supply on site

The installation of the power supply for connection of our devices must take place in accordance with international provisions and the regulations derived thereof.

Our devices are constructed according to VDE protection class I and must be connected to a socket with protective conductor contact.

Environment

The machine is fabricated from reusable or recyclable material. In accordance with the directive 2012/19 Waste Electrical and Electronic Equipment legislation (WEEE), SCHEU-DENTAL accepts to take back the device under the condition of free return delivery. Alternatively, the machine can be sent to a designated collection point for the recycling of waste electrical and electronic equipment.

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managing director of SCHEU-DENTAL GmbH