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Instructions for fabricating an aligner splint with VECTOR® 30 screws in conjunction with an archwire made of CHROMIUM wire

Compared to the classic CA® set-up technique, the treatment of malocclusions with aligner splints is even more economical, faster and more effective thanks to the use of VECTOR® 30 special screws. The individual fabrication of an archwire made of Chromium wire* fh 0.8 mm (REF 8333): The wire should be bent sensitively by hand (Fig. 5) so that the metal structure has the opportunity to stretch on the outside of the bend without tearing. You should also avoid damaging the surface by applying too much pressure to the wire with pliers. Finally, the adjusted wire ends of the bend can be deburred with a standard rubber wheel. Relevant CA® Book Volume 2 Advanced, chapter: VECTOR® 30 (English: REF 5309/German: REF 5307).

Working Material:

- REF 2484 VECTOR® 30
- REF 8333 Chromium wire
- REF 3207 ISOFOLAN®
- REF 3415 DURAN® 1.0 mm
- REF 3259 DURASPLINT® monomer
- REF 3258 DURASPLINT® polymer
- REF 5302 CA® Cyano Veneer Fast
- REF 3381 DIMO® PRO
- REF 3226 BIOPERM-Trimmer

Additionally recommended material:

- REF 3442 SIL-KITT red
- REF 5300 CA® pipettes
- REF 3220 Blocking out putty
- REF 2041 CARAT® flat pliers
- REF 5372 Set-Up manual saw
- REF 5263 CA® grinding set
- REF 3712 SD polishing brush
- REF 5288 Diamond separating disc
- REF 3714 SD high gloss polishing disc
- REF 3710 SD polishing paste A
- REF 3606 SD separating disc, serrated
- REF 3460 SD foil scissors A
- REF 3461 SD foil scissors B

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
Cut off ISOFOLAN® at model rim and insert interdental cuts using a scalpel.

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
Embed model into the pellets and thermoform DURAN® 1.0 mm. Then roughly cut it out using the SD foil scissors B.

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Push the two VECTOR® 30 screws onto the CA® arch CHROMIUM.

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
For easier insertion the ends of the CA® arch CHROMIUM are conically tapered.

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
Reduce CA® arch CHROMIUM and bend the retention loops, making sure the retention ends can be used as selective supports on the splint in the molar region.

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Separate the material within the marked area using the set-up manual saw or separating disc. If necessary, leave the splint on the model during this process for stability reasons.

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If necessary, the lower part of the plastic space holder can be reduced, in order to make sure the screws are positioned in a right angle on the anterior teeth.

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Lift off the splint from the model and combine/extend the saw cuts between the screws with the HM tricutter fine (green).

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
Fix the CA® Arch CHROMIUM and the VECTOR® 30 screws with CA® Cyano Veneer Fast on the DURAN® material. Apply just one drop each to the opening of the bent retentions using the CA® pipette (see fig.)

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Shorten the splint using SD foil scissors A and finish it with BIOPERM trimmer.

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Cover selected areas on the splint using blocking-out putty for protection.

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A dressed DIMO® PRO (see fig.) and SD high-gloss polishing mop in combination with SD polishing paste A are best suited for finishing the rims. Ensure to use them at a low speed (5.000 to 10.000 rpm).

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Finally conventionally polish the areas where you applied the resin.

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
Then roughly cut out the resin using the CA® carbide bur and mark the saw cuts with a pencil.

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Finished CA® CLEAR ALIGNER splint for the oblique segment expansion according to Dr. Pablo Echarri.

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rotation: 0.35 mm

rotation: 0.0875 mm

Note: Serious incidents in connection with the use of the product shall be reported immediately to the manufacturer and the competent authority.

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Training courses and seminars.

With the support of our internal and external speakers we offer various courses for advanced training. The training laboratory of our Academy is equipped with cutting-edge technology and can accommodate eight participants. Working in small groups guarantees intense learning combined with individual advice and support.

Our current seminars and training courses for practices and laboratories can be found at: www.scheu-academy.com

CE

SCHEU-DENTAL GmbH
www.scheu-dental.com

Am Burgberg 20
58642 Iserlohn · Germany

phone +49 2374 9288-0
fax +49 2374 9288-90

SCHEU

DE/GB 200/01/24 G REF PM0160.01