


The GEORGE GAUGE™* Bite Registration

The GEORGE GAUGE™* Bite Registration provides quick, safe and accurate data needed for the construction bite as well for intraoral anti snoring appliances as for functional orthodontic appliances.





By means of the millimeter-scale you'll have ready at hand exact and reproducible data for individual protrusion of the corresponding patient. This will serve as a base for the optimum protrusion for the construction hite


Additionnally, by means of the bite fork available in 2 and 5 mm for one-time use you can determine the necessary bite registration.

* The GEORGE GAUGE™ is patent Pending, and a Trade Name of Dr. Peter T. George, Honolulu, Hawaii, Pat ges. DE 42 92 389.1-23, © Copyright and All Rights Reserved.

Delivery programme:

// GEORGE GAUGE™ Bite Registration, comes with 3 bite forks (1 for each size)	1 Set	#5471
GEORGE GAUGE™ bite fork 2 mm, grey, small (S)	10 pcs.	#5472
GEORGE GAUGE™ bite fork 5 mm, white, small (S)	10 pcs.	#5473
// GEORGE GAUGE™ bite fork 5 mm, white, large (L)	10 pcs.	#5462







Instruction for use


- Body of GEORGE GAUGETM
- 2 Millimeter scale
- 3 Lower turn screw
- 4 Upper turn screw
- 5 Lower midline indicator
- 6 Lower incisor notch


- 7 Lower incisor clamp
- 8 Bite fork
- 9 Prongs of bite fork
- 10 Upper midline indicator
- 11 Upper incisor notch
- 12 Marking end of bite fork











- ① Loosen lower turn screw ③ and slide lower incisor clamp ⑦ forward. Center lower midline indicator ⑤ over lower central incisor, cinch up lower incisor clamp ⑦ and tighten lower turn screw ③. Remove from mouth and place bite fork ⑧ into body ①, of GEORGE GAUGE™. Use gray fork for 2 mm between incisors or white for 5 mm.
- ② Return GEORGE GAUGE™ to mouth with lower incisor notch ⑤, centered over lower incisors, and instruct patient to close into upper incisor notch ⑪, with upper midline indicator ⑩, between upper incisors. Use acrylic bur to modify upper incisor notch ⑪, if upper incisors badly rotated. Instruct patient to slide mandible first in centric relation, then into full protrusive as you observe these positions on millimeter scale ②.
- (3) From this protrusive range calculate amount of protrusion needed for the appliance you are constructing for this patient. Remove from mouth, and set marking end of bite fork (12), over appropriate position on millimeter scale (2), and tighten upper turn screw (4).
- ④ Please apply bite-taking wax on prongs of bite fork ⑤. Return GEORGE GAUGE™ to mouth with lower notch centered over lower incisors. Hand patient mirror, and instruct to close into upper incisor notch ⑥. If technique calls for maintaining a midline discrepancy, place mark on upper incisor to guide patient to proper transverse closure.
- (5) After registration material has sufficiently hardened remove from mouth. Send construction bite and bite fork (3), along with models to lab.

Further hints:



As an alternative for bite-taking wax you may as well use quick-setting, addition curing bite registration material on the base of vinylpolysiloxane (a-silicone). In this case make sure the retentions of the bite prongs are well coated.

We recommend disinfectant spray. Don't use liquids containing formaldehyde. Shouldn't be exposed to temperatures > 95°C / 200°F.

 ϵ

