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## OSAMU-Retainer® - Hints for Fabrication Procedure

In order to pressure mould BIOPLAST® 1,5 mm we recommend to use only flat trimmed models to be placed on the model platform with the even side.

For absolute transparency, we recommend the use of BIOPLAST® Insulating Agent (REF 3189).

For stability reasons, trim the pressed BIOPLAST® mould first of all to the bottom of the model.

For trimming, you can either use scissors, scalpel, diamond cutting disk or HM-cutting bur (REF 3214).

Avoid any dirt on the BIOPLAST® sheet; before laminating BIOPLAST® with IMPRELON® S pd 0,75 mm make sure to clean the contact faces which will be laminated together with lukewarm water and dry afterwards.

For pressure moulding IMPRELON® S pd 0,75 mm embed the model in the pellets in order to reduce stretching and to guarantee the necessary stability of the sheet.

Add OSAMU-Bond with a brush on the BIOPLAST® sheet during the last 10-15 seconds of heating time of IMPRELON® S pd 0,75 mm.

For long permanent compound, we recommend to prolong the cooling time of IMPRELON® S pd 0,75 mm and to cool the OSAMU-Retainer® 4-5 minutes under pressure.

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