IMPRIMO® LC Gingiva

Produktinformation

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IMPRIMO® LC Gingiva

Product information

Description:	IMPRIMO® LC Gingiva is a light-curing resin based on methacrylate and specifically engineered for the fabrication of gingival masks; photo-polymerization is done by means of image projection systems (385 nm). IMPRIMO® LC Gingiva should be used in building processes with a z-resolution of 25-100 µm. Printer parameters available for Asiga MAX™ and Asiga PRO 4K.
Processing instructions:	Make sure IMPRIMO® LC Gingiva is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Dental molded parts can be coated and repaired in the conventional way. Leave the material to stand for approx. one hour to avoid any air inclusions.
Contra-indications:	IMPRIMO® LC Gingiva should not be used for any purpose other than the generative CAD / CAM fabrication of gingival masks. IMPRIMO® LC Gingiva is a laboratory product.
Post-processing:	Upon completion of the printing process, the objects have to be removed from the model platform. Unpolymerized resin residues can be removed completely using isopropanol and a soft brush or may be resolved in the cleaning device RS wash or IMPRIMO® Clean using IMPRIMO® Cleaning Liquid.
Post curing:	Final and complete curing of the objects is done in the curing unit RS cure or IMPRIMO® Cure. The use of a vacuum or protective gas to prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. The polymerization process in the light oven RS cure in a vacuum takes 5 minutes. The polymerisation in the IMPRIMO® Cure takes 5 minutes in a nitrogen environment. The corresponding parameters are stored in the internal programme of the RS cure and the IMPRIMO® Cure unit. Any deviation from the described post-exposure process can lead to modifications of the mechanical properties or the material colour.
Shelf life / storage:	Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.
Hazard warnings:	Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Gingiva might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Gingiva. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap. Attention: Polymerized resins are chemically resistant, avoid stains on clothes.
Disposal:	According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point. All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.
Technical data:	For professional use only. Density: approx. 1.1 g/m Viscosity: approx. 2.5 Pa s Cured material: (depends on irradiation unit) Elongation at break: 90% Tear strength: 5 MPa Shore hardness: 50 A *These data result from measurements of a representative sample. Contains: urethane methacrylates, initiators, stabilizers and dyes. Material Safety Data Sheet: www.scheu-dental.com/en/download

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