UNIQ®/UNIQ®+ pd

Tiefziehfolien
thermoplastic foils
plaques de thermoformage
fogli termoformanti
planchas de termomoldeo
películas de termoformagem
giliojo tempimo plėvelės
dyptrekkingsfolier
filme de termoformare
hlbokoťažné fólie
basınçlı kalıplama folyoları
djupdragningsfilm
dybtræksfolier
syvävetokalvot
термоформовъчни фолиа
termoplastické fólie
φύλλα νάρθηκα
mélyhúzó fóliák
folie do termoformowania
termoformirajuće folije
Termoformēšanas plēves

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GA 2.000/11/26 G REF PM5385.01_Rev. 04

EN Instructions for use

Intended purpose: Thermoplastic foils for the manufacture of appliances and splints for dental and orthodontic treatments.

Indication: Suitable only for the fabrication of aligner splints (for correcting tooth and jaw misalignments) and retention splints. Instructions for use can be found at www.scheu-dental.com/pressure-moulding-material.

Contraindication: Allergic reactions to components of UNIQ®/UNIQ®+ pd (copolyester) cannot be ruled out in sensitised people.
Any use not stated above is contraindicated.

Patient target groups: Adults, Adolescents

Storage: Store at 5°C to 35°C. Protect from sunlight. Keep dry.

Note: Must only be used by trained personnel in the dental field and within 15 minutes after opening. Only open the barrier bag if it is undamaged and still sealed. For single use. Not suitable for constructing splints with acrylic resin. Any serious incident involving the product must be reported to the manufacturer and the competent authority immediately!

Safety information: The practitioner is responsible for correct use. The splint must not be used if the material is cracked/broken/contains bubbles; there is a risk of swallowing or aspiration. When using UNIQ® with 3D-printed models, a non-perforated protective foil (ISOFOLAN®) or alternatively UNIQ®+ pd must be used. Remove the protective foil after the thermoforming process. Do not bring the products into contact with alcoholic/organic solvents.

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