

Drying Oven



DURASOFT®

IMPRELON® "S"



Many acrylics are hygroscopic-reactive, i.e. they absorb humidity. This is an unpleasant property when using the pressure moulding technique as this can lead to air bubbles when plastifying materials such as DURASOFT® and IMPRELON® "S" under the infrared-heater.

In order to avoid this negative effect it is necessary to predry such high-quality materials.

The maintenance-free drying oven is perfect for it.

We guarantee its reliability.

Now this procedure can easily be done with the Drying Oven in an economic and energysaving (60 W) way.

Up to 17 plates can be predried simultaneously at a constant temperature of appr. 80°C.

Please pay attention to more de-tailed information on the added sheets in each package.

Technical Data:

Voltage: 230 V

Power: 60 W

Dim.: 24 x 20 x 11 cm

Weight: 3,5 kg