

## PRODUKT-DATENBLATT

### - Zangen CA<sup>®</sup> TIP -

CA <sup>®</sup> TIP 1	CA <sup>®</sup> TIP 2	CA <sup>®</sup> TIP 3	CA <sup>®</sup> TIP 4
Art.-Nr. 5281	Art.-Nr. 5291	Art.-Nr. 5292	Art.-Nr. 5295

Material:

Werkstoff: AISI 420 (Edelstahl)

#### Legierungsbestandteile AISI 420 (Edelstahl)

C	≥	0,15 %
Si	≤	1,00 %
Mn	≤	1,00 %
P	≤	0,040 %
S	≤	0,030 %
Cr	=	12,00 - 14,00 %
Mo	=	---
Ni	=	---
Sonstige	=	---
Fe	=	Rest