Quad Helix

Description/indication: Prefabricated expansion arch made of CHROMIUM wire (\emptyset 0,9 mm) for attachment to the upper molars in palatal sheaths.

Contraindication: Product contains nickel and chromium. In case of allergies to these components, the Quad Helix should not be used.

Application: Make sure to ligate the Quad Helix firmly into the lingual/palatal sheaths using steel or elastic ligatures.

Note: Only to be used for the above indication and by qualified personnel.

For single use. Serious incidents in connection in connection with the use of the product must be reported immediately to the manufacturer and the authority.





C€0044

SCHEU-DENTAL GmbH www.scheu-dental.com Am Burgberg 20 58642 Iserlohn · Germany phone +49 2374 9288-0 fax +49 2374 9288-90

