

# Product datasheet SALIPROTECT®

# **1** Identification of the substance / preparation / company

### **Product details**

Trade name:

Application of the substance / preparation: Manufacturer / Supplier: SALIPROTECT<sup>®</sup>

Manufacture of dental pressure moulding splints. Scheu Dental GmbH Am Burgberg 20 58642 Iserlohn Germany Fon. 0049 2374 9288-0

## 2 Composition / information on ingredients

**Chemical characterization** CAS-Number: Designation:

24937-78-8 Ethylenvinylacetat (EVA)

### 3 Handling and storage

Handling Information for safe handling:

When using do not eat, drink or smoke. Provide suction extractors if dust is formed.

### Storage

Storage conditions:

Store dry and dark at max. 20°C

## 4 Physical, chemical, mechanical and biological properties

### 4.1 General properties

Properties	Guideline	Value
Form	-	Solid
Colour	-	Transparent
Odour	-	Odourless
Density	ISO 1183	0,96 g/cm <sup>3</sup>
Water absorption after 24 h at 23°C	ISO 62 Methode 1	< 0,2 %
UV-stabilized	-	yes

Revision: 2 36/13 Replaced: 1 10/13	Page 1 von 2



# Product datasheet SALIPROTECT®

4.2 Mechanical properties				
Properties	Guideline	Value		
Tensile strength	ISO 527	>18 Mpa		
Flectional strength	ISO 527	-		
Impact strength 23°C	ISO 179/1eU	no fracture		
Notched 23°C	ISO 179/1eA	no fracture		
Yield strain	ISO 527	-		
Yield stress	ISO 527	13 MPa		
Elongation at tear	ISO 527	800-900 %		
E-modulus	ISO 527	15 MPa		
Hardness Shore A	DIN 53505	85		
Hardness Shore D	DIN 53505	-		
Rockwell Hardness	ISO 2039-1	-		

# 4.2 Mechanical properties

### 4.3 Thermal properties

Vicat softening point	ISO 306 Methode B	42°C
Thermoform resistance	ISO 75 Methode A	70°C
Continuous stress temperature		60°C

### 4.4 Biological properties / Biocompatibility

The material has been tested on biocompatibility according to DIN EN ISO 10993. It meets the requirements regarding biological compatibility for medical products / devices.

### **5 Stability and reactivity**

#### Thermal decomposition / Conditions to avoid:

No decomposition by indended use.

### 6 Disposal

The material can be recycled after separation or deposed of like commercial or household waste.

The aforementioned data are given most conscientiously but without any obligation. Any processing details are provided merely for guidance: it is the user's responsibility to check the suitability of the product for the intended application.

Revision: 2 36/13	Page 2 von 2
Replaced: 1 10/13	