

Product datasheet Lingual Bar

1 Description of the product / intended use / company

Product group:	Bars
Trade name:	Lingual Bar
Intended use:	Manufacture of removable prosthetic, orthodontic and general appliances used in the patient's mouth for the correction of tooth and jaw malocclusions or for prosthetic tooth replacement.
Manufacturer / Supplier:	SCHEU-DENTAL GmbH Am Burgberg 20 58642 Iserlohn Germany Tel. +49 2374 9288-0 gmra@scheu-dental.com

2 Composition / information on ingredients

Chemical characterization

Designation:	Stainless steel 1.4310 X10CrNi 18-8
--------------	--

Components:

Description	CAS No	Quantity
Carbon (C)	7440-44-0	= 0.05 – 0.15 %
Silicon (Si)	7440-21-3	≤ 2.00 %
Manganese (Mn)	7439-96-5	≤ 2.00 %
Phosphorus (P)	7723-14-0	≤ 0.045 %
Sulfur (S)	7704-34-9	≤ 0.015 %
Chromium (Cr)	7440-47-3	= 16.00 – 19.00 %
Molybdenum (Mo)	7439-98-7	≤ 0.80 %
Nickel (Ni)	7440-02-0	= 6.00 – 9.50 %
Vanadium (V)	7440-62-2	---
Miscellaneous	7727-37-9	N ≤ 0.11 %

Product datasheet Lingual Bar

3 Handling and storage

Handling

Information for safe handling:

Only to be used for the above-mentioned indication and by qualified specialist personnel.

Contains nickel; allergic reactions in sensitized persons cannot be ruled out

Not suitable for MRI.

Parts of orthodontic appliances can break off and be accidentally inhaled or swallowed. The appliance must no longer be worn if it is visibly damaged.

Storage

Storage conditions:

No special conditions required

4 Physical, chemical, mechanical and biological properties

4.1 General properties

Properties	Guideline	Value
Form	-	solid
Colour	-	grey
Odour	-	odourless

4.2 Biological properties / Biocompatibility

A biological evaluation was carried out for the material in accordance with DIN EN ISO 10993-1. It meets the requirements regarding biocompatibility for medical devices.

5 Stability and reactivity

Thermal decomposition / Conditions to avoid:

No decomposition by intended use.

6 Disposal

The raw material can be recycled after separation or disposed of like commercial or household waste.

The above information is provided to the best of our knowledge. In addition to the information provided, it is recommended to check the suitability of the products individually regarding the respective application and associated procedures and purposes.

Serious incidents in connection with the use of the product shall be reported immediately to the manufacturer (gmr@scheu-dental.com) and the competent authority.