



Asiga MAX™ Low Force Build Tray

The new Low Force Build Tray for the Asiga MAX™ with double Teflon foil.
Designed for soft materials.

The new Low Force Build Tray for Asiga MAX™ is specifically designed for very elastic and soft materials, such as IMPRIMO® LC Splint flex or IMPRIMO® LC Gingiva. Thanks to its double foil technique, the tray is optimally suited for 3D-printed objects with high release forces, fine details and thin wall areas, such as bruxism splints or gingival masks.

In addition, the tray provides a matter surface of 3D printed objects compared to the conventional Asiga MAX™ tray.

The Low Force Build Tray is available as 1 l tray.

Features and benefits:

- // Reduces release forces
- // Designed for extra soft materials
- // Matter surface of 3D printed objects



Delivery program:

// Asiga MAX™ Low Force Build Tray

6520.1

Do you have questions? Visit our website
www.scheu-dental.com/sales_export or get in touch with your local dealer.

DE/GB PDF 11/21 REF PM0345.01