

P

PROFESSIONAL

N O W W I T H U N I Q F O I L



NOW FOR ALL INDICATIONS:

**THE ALIGNER SYSTEM
WITH ALL-ROUND SUPPORT**

DISCOVER PROFESSIONAL. EVERYTHING AT A GLANCE.

WELCOME TO THE PROFESSIONAL WORLD 3

The aligner system for every indication.

DISCOVER THE SCOPE OF TREATMENT 5

Orthodontics can be so comprehensive - modular, tailor-made, modern.

YOUR ADVANTAGES AT A GLANCE 6

Achieve more - with less effort. For practices that think ahead.

THE MULTIPHASE TREATMENT SYSTEM 8

Success in stages - how our system takes your treatment to a new level.

OUR PREMIUM SERVICE 10

More than just support - personal service that inspires.

NOW EVEN BETTER, PROFESSIONAL WITH UNIQ 12

More comfort. More efficiency. More safety.

ADD-ONS THAT MAKE THE DIFFERENCE 13

Attachments that work - smart, discreet, effective.

MEMOTAIN® 14

precise, comfortable, and efficient.

MEMOTAIN® Peek 3D 15

3d design, metal-free, and invisible.

HOW TO ORDER ONLINE 16


Simple. Fast. Digital. Your path to the product - in just a few clicks.

THE DISCOUNT SYSTEM 18

Save smart - and enjoy attractive benefits.

PROFESSIONAL. THE ALIGNER SYSTEM FOR EVERY INDICATION.

With **PROFESSIONAL**, you can treat malocclusions efficiently and precisely - from aesthetic corrections to complex cases. Thanks to a multi-phase treatment concept and first-class support, the system offers maximum control and predictability for optimal treatment results.




PROFESSIONAL: EFFICIENCY MEETS PRECISION.

Three variants, one goal. With PROFESSIONAL, we offer you a modular aligner system that is tailored to the different orthodontic requirements in your practice - from aesthetic to complex cases. Each variant is designed to deliver **practical solutions with a high degree of predictability and clinical safety.**

SMILE

Treatment time:
up to 10 weeks

For aesthetic treatments from 5-5.




Teeth	5-5
Center line	1 mm
Closure/gap	3 mm
Inclination	10°
Angulation	10°
Rotation	10°

- Crowding
- Gaps

LITE

Treatment time:
up to 20 weeks

For corrections of moderate malocclusions in the entire dental arch from 7-7.




Teeth	7-7
Center line	2 mm
Closure/gap	6 mm
Inclination	20°
Angulation	20°
Rotation	20°

- Crowding
- Gaps
- Front tooth extraction*
- Crossbite
- Open bite
- Deep bite

FULL

Treatment time:
up to 60 weeks

For the correction of complex malocclusions in the entire dental arch 7-7, including class correction and extraction (front teeth).



Teeth	7-7
Center line	> 2 mm
Closure/gap	> 6 mm
Inclination	> 20°
Angulation	> 20°
Rotation	> 20°
Distalization	< 3 mm


- Crowding
- Gaps
- Front tooth extraction
- Crossbite
- Open bite
- Deep bite
- Class correction

DISCOVER THE SCOPE OF TREATMENT.

Orthodontics can be so comprehensive - modular, tailor-made, modern.


SMILE

10 aligners plus
5 correction aligners included




LITE

20 aligners plus
10 correction aligners included



FULL

Unlimited number of aligners until
treatment completion within 5 years



YOUR ADVANTAGES AT A GLANCE.

Achieve more - with less effort.
For practices that think ahead.

COMPETENT
ADVICE

PROVEN
MOVEMENT PRINCIPLES

OPTIMAL
THERAPY CONTROL*

PRECISE
DIGITAL PLANNING

ATTACHMENTS
& ADD-ONS

MADE IN GERMANY



* Applies to PROFESSIONAL FULL

1.

COMPETENT ADVICE

Our support team consists of experienced orthodontists and master dental technicians. We have been successfully supporting your aligner treatments for more than 10 years.

2.

PROVEN MOVEMENT PRINCIPLES

The revolutionary multi-phase treatment system eliminates unnecessary mixing movements, ensuring particularly effective aligner therapy.

3.

OPTIMAL THERAPY CONTROL

The interim impressions taken after 20 aligners with **PROFESSIONAL FULL** enable economical and optimal control and give you the opportunity to adjust the treatment at any time.

4.

PRECISE DIGITAL PLANNING

The innovative image-based medical history and precise 3D computer simulation (setup) enable exceptionally precise treatment planning.

5.

ATTACHMENTS AND ADD-ONS

Add-ons can shorten the duration of treatment. Thanks to proven attachment shapes, even complex tooth movements can be performed.

6.

MADE IN GERMANY

Like every SCHEU-DENTAL product, PROFESSIONAL is also manufactured in Germany. This guarantees short distances, direct contact persons, and high quality.

THE MULTIPHASE TREATMENT SYSTEM.

Success in stages - how our system takes your treatment to a new level.

PHASE 1: **ARCH DESIGN PHASE**

Determination of the dental arch shape. Movements: expansion or constriction, proclination, and retraction.




PHASE 2: **SPACE MANAGEMENT PHASE**

Focus on gaining space.
Movements: distalization or mesialization.
Optional: ASR (approximal enamel reduction) extractions.




PHASE 3: **ALIGNMENT PHASE**

Alignment of the teeth on the dental arch.
Movements: rotation, angulation, and torque.



PHASE 4: **OCCLUSAL FINISHING PHASE**


Adjustment of occlusion and bite.
Movements: extrusion or intrusion, angulation, and torque.
Add-on: BITE RAMP.



PHASE 5: **CORRECTION PHASE**

Performing fine corrections such as overcorrections in various orientations.






WHY SHOULD MIXING MOTIONS BE REDUCED?

Combining different movements at the same time can pose risks and cause the aligners to no longer fit properly. The efficient treatment phases of PROFESSIONAL therefore avoid unnecessary mixed movements - for optimal anchoring and better treatment results.

PROVEN MOVEMENT PRINCIPLES FOR GREATER EFFICIENCY

The idea of dividing tooth movements into different phases of movement is not new. It is therefore not surprising that in the standard work on orthodontics, "Contemporary Orthodontics" by William R. Proffit, classic bracket therapy is divided into different phases. The PROFESSIONAL system makes use of these proven principles and efficiently adapts them to aligner therapy. This reduces unnecessary mixed movements and ensures optimal anchorage. The result: more efficient results and better planning reliability.



OUR PREMIUM SERVICE.

MORE THAN SUPPORT – PERSONAL SERVICE THAT INSPIRES.

Your time is valuable – and your patients deserve the best. That's why we offer not only high-quality products, but also comprehensive service that makes all the difference: personal, efficient, and tailored to your needs. From the initial idea to final delivery, we support you with technical expertise and direct communication – for optimal results without detours.

OUR SERVICE. YOUR ADDED VALUE.



PERSONAL CONSULTATION.

We don't offer anonymous standard solutions. Our experts – master dental technicians, orthodontists, and CAD/CAM specialists – are there to assist you personally. In individually arranged appointments, we discuss your treatment goals and the appropriate setup strategy. This creates certainty, saves time, and minimizes queries during the ongoing process.

HIGHLIGHTS:

- Individual goal & setup discussions
- On-site consultation in Hilden or remote
- Personal contact persons for clinical queries
- Support from specialists in technology and practice



SMOOTH PROCESSES.

We make case planning as easy as possible for you – efficient, transparent, and with a clear focus on your practice processes.

YOUR ADVANTAGES:

- Uncomplicated ordering and case creation
- Attachments and add-ons already integrated
- Customer service with dedicated contact persons
Mon-Thu: 7:45 a.m.-4:30 p.m. | Fri: 7:45 a.m.-2:00 p.m.
Tel.: +49 2104 80041-00



YOUR CONSULTATION APPOINTMENTS.

Plan strategically. Make confident decisions. Save time.

Our services focus on personal goal and setup discussions. Here, we coordinate your individual requirements with you – precisely, practically, and efficiently. This gives you a setup that fits – and saves you valuable time during implementation.

- Goal discussion: We record your individual treatment wishes and base our planning on them.
- Setup discussion: In close consultation with you, we optimize the simulation and thereby reduce correction loops.

NOW EVEN BETTER. PROFESSIONAL WITH UNIQ.

More comfort. More efficiency.
More safety.


- Excellent aesthetics: UNIQ film impresses with its excellent appearance and allows for discreet wearing.
- Excellent material properties: The mechanical properties of the aligners are very good and provide a reliable basis for treatment.
- Optimal force transmission: The foil ensures even and targeted force transmission to support planned tooth movements.
- High abrasion resistance: UNIQ film exhibits excellent resistance to wear and tear, ensuring consistent performance throughout its entire service life.
- Effective: Treatment with the new product is effective and reliably supports the planned course of treatment.
- Patient-friendly handling: The aligners are easy to clean and maintain, ensuring high patient acceptance. No changes to the system - you can continue working as usual. Only better.

ADD-ONS THAT MAKE THE DIFFERENCE.


Attachments that work -
smart, discreet, effective.

As part of treatment with the PROFESSIONAL system, we rely on proven attachment shapes for complex tooth movements. These attachments enable efficient force transmission and optimal anchoring of the respective tooth segment.

With specially developed or utilized add-ons, SCHEU-DENTAL is setting new standards in digital aligner therapy. These additional aids, which are either already built into the aligner or need to be attached by you (such as buttons for elastic bands), make it possible to correct complex tooth misalignments more discreetly, gently, and quickly.




Ideal for anchoring, rotation, and distalization



Ideal for anchoring and extrusion or intrusion.



Treatment of deep bite indications. Anterior tooth intrusion is promoted and the posterior tooth area is dislocated.



To support class corrections, rotational angulations, mesial distalization, extrusion/intrusion, and the treatment of crossbites.

MEMOTAIN®

PRECISE. EFFICIENT. COMFORTABLE.

With **MEMOTAIN®** you can rely on maximum stability and maximum comfort. Made from Nitinol with memory effect, it is ultra-thin (0.4 mm x 0.4 mm).

YOUR ADVANTAGES AT A GLANCE.



No slipping or tilting
thanks to innovative square design



Simple application
with transfer jig can be used
in a few minutes



High patient acceptance
Ultra-thin, smooth and easy to clean



Perfect fit & optimal hold
Digitally planned and precisely
manufactured



Fast availability
order online, in the shortest
time in your practice

+ Our added security: 5-year guarantee on material and design.

Order now & experience
retention in a new way!



MEMOTAIN® Peek 3D

3D DESIGN. METAL-FREE.

INVISIBLE.

A radiant smile stands for self-confidence, health and success. The **MEMOTAIN® Peek 3D**, made of PEEK (polyether ether ketone), a biocompatible high-performance plastic, sets new standards with its dimension of 0.45 mm x 0.9 mm in comfort, aesthetics and durability.

YOUR ADVANTAGES AT A GLANCE.



Almost invisible

Thanks to its tooth-colored, ivory-like appearance, the retainer is hardly recognizable



Metal-free & biocompatible

Allergy-friendly, no metallic taste - Ideal for sensitive patients



Perfect fit thanks to 3D technology

Precise, CAD /CAM-manufactured for maximum wearing comfort



Comfortable to wear

No sharp edges, gentle adaptation adaptation to the row of teeth



3D Design

Optimum fit for tight spaces space conditions and special tooth shapes




Our added security: 2-year guarantee on material and design.

Order now & experience retention in a new way!



HOW TO ORDER ONLINE.

Simple. Fast. Digital.
Your path to the product -
in just a few clicks.




YOUR DIGITAL ASSISTANT FOR PATIENT MANAGEMENT.

The SCHEU-DENTAL Treatment Manager is your digital assistant that makes it easier for you to manage your patient cases and optimizes the ordering process for tailor-made products. Benefit from numerous practical functions and intuitive operation.

YOUR ADVANTAGES AT A GLANCE

- **Easy ordering:** Order all SCHEU-DENTAL custom-made products directly via the Treatment Manager.
- **Comprehensive archive:** Keep track of all current and completed patient cases.
- **Online check:** Check and approve set-up forecasts conveniently online.
- **Status tracking:** Track the current status of your treatments in real time.
- **Simple data upload:** Upload 3D scan data easily and manage it centrally.



REGISTER NOW AND GET STARTED:
ORDER.SCHEU-DENTAL.COM

SHIPPING COSTS AT A GLANCE.

We rely on professional shipping to ensure that your order arrives quickly and reliably. The shipping costs are transparent and fair - you can find all the details here:

	Germany	Austria	Switzerland
Standard shipping	7,90 €	8,90 €	17,90 €
Express shipping	12,90 €	17,90 €	26,90 €

Shipping costs for other countries on request!

TERMS AT A GLANCE.

Clear timelines for maximum planning reliability: from data analysis to completion - that's how long it takes to produce your SCHEU-DENTAL custom-made products.

Production/delivery times (in working days)		Creation of the setup forecast (from receipt of the data/models)		Production (from release of the setup forecast)	Shipping
		Data	Models		
PROFESSIONAL Smile	Aligner STL data	4 days	6 days	1 day	-
	Aligner splints	4 days	6 days	5 days	approx. 2 days
PROFESSIONAL Lite	Aligner STL data	5 days	7 days	1 day	-
	Aligner splints	5 days	7 days	5 days	approx. 2 days
PROFESSIONAL Full	Aligner STL data	6 days	8 days	1 day	-
	Aligner splints	6 days	8 days	5 days	approx. 2 days

The processing times for subsequent steps and correction packages depend on the package selected. If a new impression or scan is required, the above planning and production times apply.

ORDER. POINTS. PROFIT.

YOUR ADVANTAGE: SAVE UP TO 47 %!

Take advantage of our multi-level discount system and secure attractive discounts! The more treatments you carry out with SCHEU-DENTAL custom-made products, the greater your savings.

- **Six discount levels:**
depending on the points you have collected.
- **Annual revaluation:**
Your discount level is determined each year on 30.06. based on the previous 12 months. On 01.07. your points counter is set to zero.
- **Flexible upgrade:**
If you reach a higher level in the current year, your discount will be adjusted on a weekly basis.
- **Discount on all treatment products:**
Benefit automatically from better conditions!

OUR DISCOUNT LEVELS FOR YOU!

Discount level point	I ≤ 300	II + 300	III + 800	IV + 2.000	V + 5.000	VI + 15.000
Aligner	0 %	5 %	10 %	20 %	32 %	47 %
Retainer	0 %	5 %	10 %	15 %	18 %	20 %

TERMS AND CONDITIONS FOR THE SCHEU-DENTAL custom-made PARTNER PROGRAM

The following conditions apply to the acquisition of discounts, the redemption of discount vouchers, and the provision of other services:

1. Customers of SCHEU-DENTAL custom-made GmbH who have a customer account in Treatment Manager are eligible to participate. There is no legal right to participation. SCHEU-DENTAL custom-made reserves the right to refuse participation in the discount program if there are legitimate reasons (e.g., abuse or late payment).
2. Participation begins after the cut-off date on July 1, 2021. The customer's points from the last twelve months are calculated on the cut-off date (June 30). The customer is then assigned to the respective discount level. July 1 is the start date. From this date, all customers will receive at least the guaranteed discount according to their discount level, and the points counter will be reset to 0. Customers who open a new customer account after the cut-off date will start at the lowest discount level with a counter of 0.
3. The discount program is based on the total number of points assigned to specific products when ordering custom-made treatments from SCHEU-DENTAL. The amount of the discount credit is divided into six levels and is granted when a discount level is reached.
4. Discounts in accordance with the applicable discount level are granted on all treatment sales. This does not include shipping costs, cancellation fees, other non-discountable sales, and individual products that are excluded from the discount.
5. Discount levels are updated at the end of each month. At the end of each month, a server query is performed on the customer account. The new discount rate is then applicable to the customer on the following day.
6. Cannot be combined with other promotional discounts.
7. Conversion into cash is not possible.
8. Customers can check their current discount level by calling customer service during service hours. SCHEU-DENTAL custom-made strives to ensure the highest possible level of data security. Nevertheless, SCHEU-DENTAL custom-made is not liable for damage to customer data caused by malware or errors in the software.
9. The discount level cannot be transferred to third parties. Indirect transfer of the discount level by placing orders in one's own name for third parties is prohibited. In the event of misuse for which the customer is responsible, SCHEU-DENTAL custom-made reserves the right to refuse to honor the discounts and/or to exclude the customer from participation in the discount program.
10. SCHEU-DENTAL custom-made reserves the right to discontinue, modify, or replace the discount program with another program after giving reasonable notice, or without notice if there is an important reason to do so. Points or discount entitlements already acquired by participants up to that point will remain unaffected and can be redeemed within a reasonable period of time in accordance with the previous conditions.
11. Personal data will only be processed and used for purposes related to the implementation of the program. Personal data will not be disclosed to external third parties unless SCHEU-DENTAL custom-made is required to do so by law or court order. Otherwise, personal data will be deleted after the expiry of the commercial and tax law retention periods, unless there are legal reasons to the contrary. With regard to the stored personal data, the customer is entitled to the statutory rights of information, correction, blocking, and, if necessary, deletion.
12. Reference is made to the law on combating corruption in the healthcare sector.
13. Should one or more provisions of this contract be or become legally invalid, this shall not affect the validity of the remaining provisions. The same applies in the event of a loophole. These terms and conditions correspond to the status as of August 19, 2021.

If you have any difficulties understanding this information, SCHEU-DENTAL custom-made customer service will be happy to help you on +49 2104 80041-00.

Our general terms and conditions apply. These can be viewed at <https://www.scheu-dental.com/agb>



ARE YOU EXCITED TOO?
DISCOVER MORE. SCAN NOW.

YOUR DIRECT LINE TO US!

Do you have questions about your order, shipping or would you like individual advice? Our team is here for you!



**Info-Hotline /
Distribution**

02104 800041-00

info@scheu-dental.com



**Clinical support
Retainer**

02104 800041-01

memotain-support@scheu-dental.com



**Clinical support
Aligner**

02104 800041-02

aligner-support@scheu-dental.com

OUR SERVICE TIMES

Monday - Thursday 07:45 - 16:30
Friday 07:45 - 14:00



SCHEU-DENTAL custom-made GmbH
Walder Straße 53, 40724 Hilden, Germany

phone +49 2104 80041-00
fax +49 2104 80041-99

E-Mail info@scheu-dental.com
Web www.scheu-dental.com