

dreamTAP[™]-splint

An innovative step forward in the TAP[®] system

The treatment of snoring and sleep apnoea with intraoral appliances such as the TAP[®] splint has proven successful as independent therapy in sleep medicine as well as supplemental treatment to respiratory masks (nCPAP). Highest wearing comfort and reliable functionality of the TAP[®] splint are guaranteed by its constructive advantages, such as continuous adjustment and lateral excursion.



- // nickel-free chromium-cobalt alloy
- // custom-made
- // infinitely variable
- // 15 mm of protrusion range thanks to exchangeable hooks
- // 10 mm of lateral excursion
- // extremely break-resistant
- // outstanding wearing comfort
- // long-standing clinical experience

Based on proven experience: The new dreamTAP[™] device

The dreamTAP[™] device offers even greater wearing comfort thanks to the modified insertion position. Additional tongue space is created by the protrusion adjustment device in the lower jaw. With its exchangeable hooks, the dreamTAP[™] device allows a protrusion range of up to 15 mm.

The dreamTAP[™] components are metal injection moulded using a nickel free chromium-cobalt alloy and ensuring outstanding stability.



Manufactured by
AM AIRWAY
MANAGEMENT

CE

DE/GB 300/11/24 G REF PM0193.01

