

REF 3214 - Cutting bur (white) Rough cutting of hard foils	max. 25.000 r/min
REF 3369 - HM-carbide cutter (yellow) Rough cutting of hard foils and hard/soft foils	max. 7.500 r/min
REF 3370 - HM-Tricutter fine (green) Finishing of hard foils and hard/soft foils	max. 7.500 r/min
REF 3377 - Finishing bur (red) Cutting of hard foils and splints with integrated wire elements	max. 7.500 r/min
REF 3226 - Bioperm-Trimmer (blue) Rough trimming of soft foils and hard/soft foils	max. 25.000 r/min
REF 3371 - Polyfix (brown) Finishing of soft foils and hard/soft foils	max. 12.000 r/min
REF 3247 - OSAMU-Polisher (grey) Final polishing of hard/soft foils	max. 15.000 r/min
REF 3376 - DIMO®PRO SLIM (orange) Final trimming and polishing	max. 10.000 r/min

DE/GB 500/03/25 G REF 0306.01