CARAT® Wire Bending Pliers

Pliers made of high-grade stainless steel alloy with exactly machined tips (depending on type from 1.0 to 1.5 mm) and precisely guided pivots. Ergonomic design with rounded edges ensures fatigue-free handling.

CARAT® Krampon Pliers

Pin and wire cutting/bending pliers with chequered and crenate tips, for hard wires up to \emptyset 1.2 mm and soft wires up to \emptyset 1.5 mm.

○ CARAT[®] Krampon Pliers

#2040

CARAT® Flat Pliers

Flat needle-nose pliers acc. to Marburg model with delicate, chequered guiding tips.

○ CARAT[®] Flat Pliers

#2041

CARAT® Young Pliers

Loop and arch bending pliers with multiple creneate tips for labial arches and hard wires (\emptyset to 0.7 mm).

○ CARAT[®] Young Pliers #2042

CARAT[®] Angle Pliers (HM)

Carbide wire bending pliers with a square end and a round tip, for wires up to 0.7 mm.

○ CARAT[®] Angle Pliers (HM)

#2043

CARAT® Aderer Pliers large

Aderer wire bending pliers, three prong with reinforced tips, for wires up to 0.9 mm.

○ CARAT[®] Aderer Pliers large

#2044

CARAT® Aderer Pliers medium

Aderer wire bending pliers, three prong with reinforced tips, for wires up to 0.7 mm.

○ CARAT[®] Aderer Pliers medium

#2046

phone +49 23 74 92 88-0 fax +49 23 74 92 88-90 Am Burgberg 20 58642 Iserlohn · Germany

ermany

S

Dental Technology

1.0 mm

1.0 mm

26.5 mm







1.5 mm

C