

# Stainless Steel Holder

## Instructions for use

**Description:**

Holder made of stainless steel (Mat. No. 1.4305 / AISI 303) with straight and angled holder.

**Indication/****contraindication:**

For safe use of CA® Stripping Tools without EVA handpiece. Also suitable for all instruments with a shaft diameter of 2.35 mm.


**Intolerances:**

Contains nickel; may cause allergies in sensitised individuals.

**Attention:**

To ensure a secure hold of the instruments used in the stainless steel holder, ensure to tighten the union nut using suitable pliers, if necessary.

To be used only for the above indication and by qualified personnel. Serious incidents in connection with the use of the product must be reported immediately to the manufacturer and the competent authority!



Sterilizable

## REF 5358



SCHEU-DENTAL GmbH  
www.scheu-dental.com

Am Burgberg 20  
58642 Iserlohn · Germany

phone +49 2374 9288-0  
fax +49 2374 9288-90

**SCHEU**

DE/GB 500/01/24 G REF PM0253.01

