

Aligner Pliers (CA® TIP 1-4)

Description:

CA® TIP 1

Pliers for aligners with large spherical tip allow insertion of retention points in CA® splints, for increase of friction.

CA® TIP 2

Pliers for aligners with small spherical tip allow forming of activation points in CA® splints.

CA® TIP 3

Pliers for aligners with oval plateau, enable creation of activation surfaces in CA® splints.

CA® TIP 4

Pliers for aligners destined for torque control with curved plateau, enable creation of labial / lingual activation surfaces in CA® splints.

Cleaning/Care:

We recommend thermal disinfection (delicate program not exceeding 121°C) or spray disinfection.

Warranty:

5 years on material, workmanship, breakage and corrosion.
All pliers are made of stainless steel.

Hint:

To be used only for the above indication and by qualified personnel. For single use. Serious incidents in connection with the use of the product must be reported immediately to the manufacturer and the competent authority!

Delivery program:

CA® TIP Pliers		
CA® TIP 1	5281.1	1 piece
CA® TIP 2	5291.1	1 piece
CA® TIP 3	5292.1	1 piece
CA® TIP 4	5295.1	1 piece

Detailed information on use of the pliers can be found here:



SCHEU-DENTAL GmbH
www.scheu-dental.com

Am Burgberg 20
58642 Iserlohn · Germany

phone +49 2374 9288-0
fax +49 2374 9288-90

SCHEU
SCHEU GROUP

