

IMPRELON® S pd / IMPRELON® S+ pd

Now
also available
with protective
film!



The ideal solution for invisible splint retainer.

Our IMPRELON® S is the ideal solution to be used as a retention foil. The improved material properties convince with extremely break- and abrasion-resistancy. The foils with a thickness of 0.75 or 1.0 mm and a diameter of 125 mm are pre-dried and individually vacuum-packed for immediate use.

This foil is now also available as a '+' version and offers additional convenience in production thanks to its integrated protective foil.



Stability

IMPRELON® S provides a constant force with improved retention by holding the teeth longer in position.



Resistance

The improved material of IMPRELON® S pd is more resistant to cracks and scratches.



Comfortable

IMPRELON® S ensures a high and constant restoring force during continuous stress and relief, promoting consistent tooth movement and enhanced wearing comfort for the patient.




Variety

Whether round or strip, standard or '+' version, IMPRELON® S is available in numerous designs, offering maximum flexibility.

High and constant force with extra comfort.


IMPRELON® S provides a constant force with improved retention by holding the teeth longer in position. Here in comparison with PETG and a competitor's material.

Constant spring force for improved retention



Stability that lasts:

IMPRELON® S impresses in endurance testing



Dimensions

IMPRELON® S pd	125 mm Ø	Pack size
0.5 mm	REF 3320	10 pcs
0.75 mm	REF 3322	10 / 100 pcs
1.0 mm	REF 3324	10 / 100 pcs
1.5 mm	REF 3326	10 pcs
2.0 mm	REF 3328	10 pcs
3.0 mm	REF 3330	10 pcs

CE 0044



IMPRELON® S+ pd	125 mm Ø	499 x 114 mm	Pack size
0.625 mm		REF 3334	10 / 100 pcs
0.75 mm	REF 3331	REF 3335	10 / 100 pcs
1.0 mm	REF 3332	REF 3336	10 / 100 pcs

CE 0044



Do you have any questions?

Contact our technical customer service team or make an appointment with our field service: <https://www.scheu-dental.com/en/contact>

SCHEU-DENTAL GmbH
scheu-dental.com

phone +49 2374 9288-720
fax +49 2374 9288-90

