



ANNIVERSARY PROMOTION

UNIQ pd/UNIQ+ pd

10 + 2 = 100*

NOT EVERY CALCULATION ADDS UP.

**THIS ONE DOES:
A CAMPAIGN TO CELEBRATE 100 YEARS OF SCHEU.**

*As part of our anniversary promotion, when you purchase 10 packs of UNIQ pd/UNIQ+ pd, you will receive two additional packs free of charge. Promotion period: April to June 2026

ANNIVERSARY PROMOTION

BUY 10 PACKS.
GET 12.



10+2

UNIQ pd/UNIQ+ pd

THE IDEAL SOLUTION FOR EVERY ALIGNER TREATMENT.

Our innovative thermoforming foil sets new standards for greater efficiency and comfort in aligner treatment. UNIQ pd/UNIQ+ pd achieves with just one layer what many multi-layer foils cannot. The single-layer material enables it to withstand forces throughout the entire treatment. The result? Comfortable to wear with gentle initial forces and constant restoring forces. With UNIQ pd/UNIQ+ pd, we have developed a foil that achieves maximum benefits with just one layer for every aligner treatment!



One layer, maximum effect

Thanks to its improved material properties compared to conventional single-layer aligner materials, UNIQ pd/UNIQ+ pd exerts less pressure while maintaining it for a longer period of time. For the patient, this means gentler and more effective tooth movement combined with greater wearing comfort.

General information	UNIQ pd/UNIQ+ pd	
Thicknesses (in mm)	UNIQ pd:	UNIQ+ pd:
0.5 mm	3647	3651
0.625 mm	3648	3652
0.75 mm	3649	3653
Recommended foil thicknesses for aligner therapy 0.625 mm.		

For insulation, we recommend our ISOFOLAN® foil. The '+' version is already equipped with a workable insulation foil.

Integrated protective foil

Our UNIQ pd/UNIQ+ pd foil has a protective foil already integrated. This saves you time and material when thermoforming without compromising on precision.



Safely packaged and ready for immediate use

Safely packaged and ready for immediate use
UNIQ pd/UNIQ+ pd foils are pre-dried and individually vacuum-packed. This provides with special protection against external influences and extends the storage life.

Discover our Aligner Family!

As part of the Aligner Family, UNIQ pd/UNIQ+ pd foil is ideally suited for use in aligner therapy. Discover other foils from our range and find the right foil for every application with SCHEU-DENTAL.

